IKA® MF 10

Betriebsanleitung	DE	5
Operating instructions	EN	16
Mode d'emploi	FR	27
Instrucciones de manejo	ES	38
Handleiding	NL	49
Istruzioni per l'uso	IT	60
Driftsanvisning	SV	71
Driftsinstruks	DA	82
Driftsvejledning	NO	93
Käyttöohje	FI	104
Instruções de serviço	PT	115

Υπδείξεις ασφάλειας	EL	126
Wskazówki bezpieczeństwa	PL	127
Bezpečnostní upozornění	CS	128
Biztonsági utasítások	HU	129
Varnostna navodila	SL	130
Bezpečnostné pokyny	SK	131
Ohutusjuhised	ET	132
Drošības Norādes	LV	133
Saugos Reikalavimai	LT	134
Инструкции за безопасност	BG	135
Indicații de siguranță	RO	136

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Contents

	Page
CE-Declaration of conformity	16
Explication of warning symbols	16
Safety instructions	17
Unpacking	18
Useful information	18
Correct use	19
Commissioning	20
Eliminating malfunctions	24
Maintenance and cleaning	25
Accessories	25
Technical data	26
Warranty	26

CE - Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2006/95/EEC and 2004/108/EEC and conforms with the standards or standardized documents: DIN EN IEC 61 010-1; DIN EN ISO 12 100-1, -2; EN 60 204-1 und DIN EN IEC 61326-1.

Explication of warning symbols

General hazard

Note!

This symbol identifies information **that is of importance for the technically correct functioning of the system**. Disregarding this information can result in damage to the appliance or to system components.

Attention!

This symbol identifies information **that is of vital importance for safeguarding your health and safety**. Disregarding this information can lead to health impairment and injuries.

Safety instructions

For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of:
 - projectile parts.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The mill is to be operated only under the supervision of an authorized operator..
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- The mill is to be switched off completely using the rocker switch (Pos. A 13) before opening the door.
- Excess pressure develops with closed collecting vessels, in particular when using dry ice or liquid nitrogen. This can lead to ground material escaping upwards out of the funnel.
- Feeding stock residue that remain in the funnel may not be manually removed during operation or with unsuitable tools e.g. screwdrivers. They should only be removed with suitable cleaning tools when the device is switched off.
- Fingers and/or other articles should not be placed into the discharge cone (Pos. P 1/ S 1) in order to clean it during operation. Cleaning is only permitted when the device is switched off and the power supply plug has been disconnected.
- No highly combustible materials may be cut up.
- **Note:** After longer grinding cycles the grinding chamber and the crushing tools may be hot. Let the mill cool off before cleaning or performing any repair work.
- **Note:** Exercise extreme caution when assembling, disassembling or cleaning the cutting-grinding head, particularly with regard to the sharp blade edges on the rotor (Pos. S 20). The-

se parts are to be handled with extreme care. Wear protective personal equipment.

- Please observe the relevant safety instructions and guidelines regarding dust explosions.
- Always wear suitable protective equipment when embrittling the mill feed with coolants (liquid nitrogen or dry ice), (see laboratory guidelines and accident prevention regulations applicable to the laboratory).
- Ensure that pressure does not build up in the grinding chamber (particularly when using coolants).
- Caution: When using liquid nitrogen, the milling process must only be started when the liquid nitrogen has fully evaporated.
- Caution: Appliance parts may still be extremely cold for some time after cooling with liquid nitrogen and should only be handled whilst wearing protective equipment.
- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA®** if you have any questions.
- Do not operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through-processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Always disconnect the plug before fitting accessories.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- After an interruption to the power supply or a mechanical interruption, the device does not start up again automatically.

For protection of the equipment

- The appliance may only be opened by experts.
- The voltage stated on the type plate must correspond to the mains voltage.
- Protect the appliance and accessories from bumps and impacts.

Unpacking

• Unpacking

- Please unpack the device carefully

- In the case of any damage a detailed report must be set immediately (post, rail or forwarder)

• Delivery scope IKA® MF 10

The delivery quantity for the drive includes:

an IKA® MF 10 drive unit, a collecting pan, a mains cable, a face spanner, an Allen wrench, a fork wrench and an operator's manual.

The delivery quantity for the impact grinding head includes:

an impact grinding surface with built-on door and three screwed-in fastening pins, a dosing sluice, a hammer rotor, a countersunk screw, a disk, a cover and the filter which you ordered.

The delivery quantity for the cutting-grinding head includes:

a cuttinggrinding head with three built-in bed knives, built-on door, tabular slide valve and funnel, as well as three screwed-in attaching pins, a plunger, a blade rotor with three mounted rotor blades, a cheese head screw, a disk and the filter which you ordered.

Useful information

The IKA® MF 10 analytical mill is a continuously operating universal mill. It makes it possible to use two different grinding methods which cover the greatest range of all grinding tasks occurring in the laboratory.

Impact grinding:

Application for hard, brittle and/or dried materials or ground materials embrittled as a result of cold (e.g. grain, coffee, coal, etc.).

In this case, the ground material is cut up by means of rotary beaters and comminuted on the grinding surface attached to the blow bars. The ground material is broken and ejected into a suitable collecting pan by means of corresponding filters. The feed dosage provides for optimal dosage so that the mill can work at the optimal operating point. The grinding chamber is locked so that it is dustproof, the

collecting pans are attached with standard ground joints. Thus the loss ground material is kept to a minimum.

Cutting-Grinding:

Application for soft fibrous ground material (e.g. hay, straw, plastic waste, paper, wood, etc.).

In this case, the ground material is cut between rotary blades and standing blades attached at the grinding surface. The ground material is conveyed into the collecting pan by means of filters here as well. The feed unit (guide channel) makes it possible to load awkward, bulky material. The cutting loads mean that less energy is introduced into the ground material and thus results in less heat generation. This keeps ground material loss to a minimum.

• Range of use

- Laboratories
- Technical colleges
- Pharmacies

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance or the printed circuit board are modified by third parties.

This device is suitable for use in all areas, except industrial areas.

• Use

The area of application for the universal mill is very extensive, the list of materials which are capable of being ground and indicated here is not exhaustive. The scope of application can be extended further by grinding attempts with the respective grinding method and setting the appropriate grinding parameters (rotor speed, filter hole diameter, sample preparation e.g. cooling).

Impact grinding:

The hammer mill crushes soft, medium-hard and brittle materials to approximately Mohs hardness 6 (feldspar: 6 Mohs, quartz: 7 Mohs). With greater hardness the wear on the grinding tools is much higher and grinding is therefore uneconomical.

Anything can be ground that breaks, is dry does not contain very fatty materials.

The following is a selection of substances drying to grind: grain, barley, corn, malt, pectin, roasted coffee, nutshells, bones, ergot, coal, charcoal, coke, peat, ash, cellulose, drugs, artificial fertilizers, refractory clay, feed, spices, glass, resin, potash, cores, salts, cinders, stones and soil (bau-

xite, gypsum, limestone, dry clay, mica, graphite, pumice stone, gneiss, magnetite, heavy spar, sandstone, brown iron ore, etc.).

The scope of application can be further enhanced through embrittlement of the ground material (e.g. with dry ice and/or liquid nitrogen). Tough ground material must be cooled, e.g. crushed dry ice put into the grinding funnel.

The ground material must be coarse-crushed into maximum 10 mm chunks beforehand.

Note: With abrasive materials (e.g. stones or glass) it will be necessary to clean the dosing sluice repeatedly because settling grinding dust in the guides and it thus rotating the dosing sluice is no longer possible.

Cutting-Grinding:

The cutting-grinding mill cuts up bulky, elastic, fibrous, cellulose and soft materials. Mixing materials such as various types of waste must be free of iron and non-ferrous metals. The feedings may not be damp and greasy, as they will otherwise stick to the grinding surface and filter.

The following is a list of the materials that may be processed: Leaves, fibers, spices, grain, hop, wood, pasteboard, paper, hay, straw, plastics, tobacco, peat, feed, pasta, roots and branches, etc.. The scope of application can be further enhanced through embrittlement of the ground material (e.g. with dry ice and/or liquid nitrogen).

Bulky material is metered with the tubular slide valve, the feed grain size should not be not greater than 15 mm. Bulky materials (e.g. plastic bags) are metered by means of plungers, but the pieces should not be too large (e.g. pieces of fir wood: max. diameter 10 mm, max. length 40 mm; pieces of foil: max. width 30 mm, max. length 100 mm).

Note!

Make certain that the mill drive stands firmly on a table. Comply with the ambient conditions given in the technical data.

The collecting pan (Pos. A 19) is placed on its base in such a way that the round magnets fit into the recess on the base.

Impact grinding mill:

Assembling the impact grinding head:

The ON/OFF switch Pos. A 13 must be set to zero and the power supply plugs must have been disconnected out before beginning assembly of the impact grinding head. First the three threaded bolts Pos. P 5 (see fig. 1) are screwed out to the same height (distance between hexagon of the bolt and the flat surface on the grinding head is approximately 2mm, this corresponds to two turns). Make sure that the flange piece on the drive and on the head are clean and free of dust. Abrading both fitting parts with a slightly oiled rag facilitates assembly. The adjustment spring must be in the area of the opening (in the back wall of the mill chamber). The head is now pushed onto the drive, while making certain that the parts do not jam. The threaded bolt ends must fit into the drillings of the transmission cover. Thus the head has been rotated slightly in a clockwise direction. After the bolts are fitted in the head is turned until it is stopped by the bolts and so that the head now stands upright. Now close and lock the grinding chamber door with the operating knob Pos. P10 (3 to 5 revolutions are sufficient). This ensures that the grinding chamber is optimally aligned so that the door can be easily opened and closed later. Now the screw mountings are tightened with the fork wrench provided. Make certain that the pins are tightened toward the head.

Note: Tighten the pins evenly (see fig. 3.1)!

Now the hammer rotor Pos. P 13 can be installed. The rotor is pushed onto the drive shaft as shown in Fig. 3, fastened with the disk Pos. P 11 and the countersunk screw Pos. P 15. The countersunk screw, P 15, is tightened with the Allen key supplied, while the rotor is held with the face spanner that is included in delivery. The spanner is positioned so that the pins of the spanner are

inserted into the bore holes on the face of the rotor. The disc must be turned so that the bore hole in the rotor is accessible for the spanner pins.

Note: After tightening the cheese head screw with the Allen key, turn the hammer rotor once (360°). The rotor may not touch. Repeat assembly if necessary.

Now the filter Pos. Z 1 that is suitable for the grinding procedure can be inserted; to do so the filter is pushed into the grinding chamber from the front. The filter and the filter holder surface in the grinding chamber must be free from dirt and dust. Now the grinding chamber door can be closed and locked with the operating knob Pos. P 10. The screw mounting of the operating knob is screwed in by pressing and twisting it into the internal thread on the drive cover. The first turns are easy, if turning becomes difficult then the door may be closed, but the safety switch is not yet unlocked. Continue turning until stop (approximately two revolutions), only then is the safety switch unlocked and the mill ready to start.

Grinding operation with impact grinding mill:

The suitable filter is pushed in and the door is locked before beginning the grinding procedure. Afterwards a collecting pan is fastened to the discharge cone (standard ground joint NS 29). Suitable containers are e.g. Erlenmeyer flasks Pos. Z 3, round-bottomed flasks, test tubes made of either glass or plastic. The containers are to be secured with a fixing clip Pos. Z 2 (see fig. 1). For heat sensitive materials or with ground materials that have a very high increase in volume during cut up, it is advisable to use a twonecked flask with a filter attached to its second (e.g. cotton wool, etc...). The ground material is cooled from the air flow and no excess pressure is created. For very critical ground materials (very fine powder after being cut up) the ground material can also be sucked off by means of a commercial vacuum cleaner with a suction tube adapted at the discharge cone. The powder is then in the vacuum cleaner bag where it can be removed for further processing steps.

The ground material which has prepared to the correct feed grain size is filled into the funnel and then the cover Pos. P 17 is applied. The speed knob Pos. A 12 is turned to the left stop (3000 rpm).

The mill is started by shifting the rocker switch Pos. A 13 to Position 1. First the green „Power“ LED is lit and then „Overload“ LED lights up briefly before the motor smoothly starts. The speed can be set by turning the speed knob to the desired speed. A speed of more than 4500 rpm is recommended for impact grinding.

Now the ground material can be metered by turning the dosing sluice Pos. P6 in a clockwise direction. You can hear if ground materials gets into the grinding chamber and if too much materials has been metered then the rotor speed decreases. This drop in speed is clearly indicated by the lower grinding noises and the orange „Overload“ LED lights up. Metering of the ground material should only be resumed when the speed has returned to the adjusted value or if the „Overload“-LED is no longer lit.

If this is done then the mill becomes overloaded and switches off after some time because of overheating. This is then indicated by the red „Error“ LED. In addition, a beep sounds to indicate that the mill should be switched off using the rocker switch Pos. A 13. The mill must then cool off for a period in order to be able to start in again (see chapter "Eliminating malfunctions").

In addition, the ground material is no longer cut up the lower speed when the unit is operated with an overload because the impact energy no longer suffices for breaking the ground material. Keeping the ground material in the grinding chamber and the ground material to heat up. In extreme cases this can result in the ground material caking on the grinding chamber parts which then entails a time-consuming cleaning operation.

If the funnel is empty then the mill remains turned until all of the ground material from the grinding chamber is in the container. This is the case if cutting up noises can no longer be heard.

When it comes to cleaning the grinding chamber experience has shown that it makes sense to briefly switch the unit to maximum speed after the grinding operation has been completed and only then to turn the speed knob to the left stop. The mill is then shut off with the switch.

The filter is removed after opening the grinding chamber door and any remaining ground material is wiped into the collecting pan.

Caution: After longer grinding cycles the grinding chamber and the filters may have an increased temperature. Protective personal equipment should be worn in this case.

It may help to crush the ground material in several steps, i.e. a filter with large perforations is used for pre-crushing and finer ones in each of the following steps. (see fig. 4).

Advantages:

- Fewer failures, e.g. due to clogged filters or overloading the motor.
- Less heat-up of ground material since the retention time of the ground material in the grinding chamber is substantially reduced.
- Closer grain size spectrum for the final material.
- Normally less total grinding with better results as despite several runs comminution is performed more quickly and more effectively.

Disassembling the impact grinding head:

Before disassembling the head of the mill, the mill itself must first cool down as the grinding chamber and the beater rotor cannot be detached otherwise.

The power supply plug must be disconnected before disassembly. After opening the door, the countersunk screw, P 15, is loosened using an Allen key, while the rotor is held by the face spanner that is included in delivery. The spanner is positioned so that the pins of the spanner are inserted into the bore holes on the face of the rotor. The disc must be turned so that the bore hole in the rotor is accessible for the spanner pins.

If the hexagon socket is dirty then it must be cleaned by means of a brush. After the countersunk screw and the disc, P 11, have been removed the hammer rotor, P 13, is removed. The adjustment spring must be in the area of the opening (in the back wall of the mill chamber). Now the three screw mountings Pos. P5 are loosened by means of a fork wrench and the grinding chamber is rotated in a clockwise direction and taken off.

Cutting mill:

Assembling the cutting mill head:

Caution: The blades on the rotor and in the grinding chamber have sharp edges. Both parts should only be handled when wearing protective gloves.

The ON/OFF switch Pos. A 13 must be to zero and the power supply plugs must have been disconnected out before beginning assembly of the cutting-grinding head.

First the three threaded bolts Pos. S 12 (see fig. 2) are screwed out to the same height (distance between hexagon of the bolt and the flat surface on the grinding head is approximately 2 mm, this corresponds to two turns). Make sure that the flange piece on the drive and on the head are clean and free of dust. Abrading both fitting parts with a slightly oiled rag facilitates assembly. The head is now pushed onto the drive, while making certain that the parts do not jam. The threaded bolt ends must fit into the drillings of the transmission cover. Thus the head has been rotated slightly in a clockwise direction. After the bolts are fitted in the head is turned until it is stopped by the bolts and so that the head now stands upright. Now close and lock the grinding chamber door with the operating knob Pos. S 9 (3 to 5 revolutions are sufficient). This ensures that the grinding chamber is optimally aligned so that the door can be easily opened and closed later.

Now the screw mountings are tightened with the fork wrench provided. Make certain that the pins are tightened toward the head (see fig. 3.2).

Note: Tighten the pins evenly!

Now the hammer rotor Pos. S 20 can be installed. For this the fitting key Pos. A 7 is first inserted into the groove of the drive shaft Pos. G 6, make certain that the groove is free from dirt and dust here as well. The rotor is pushed onto the drive shaft as shown in Fig. 2, fastened with the disk Pos. S 22 and the cheese head screw Pos. S 21. The countersunk screw, P 15, is tightened with the Allen key supplied, while the rotor is held with the face spanner that is included in delivery. The spanner is positioned so that the pins of the spanner are inserted into the bore holes on the face of the rotor. The disc must be turned so that the bore hole in the rotor is accessible for the spanner pins.

Note: After tightening the cheese head screw with the Allen key, turn the blade rotor once (360°). The rotor may not touch. Repeat assembly if necessary.

Now the filter Pos. Z 1 that is suitable for the grinding procedure can be inserted; to do so the filter is pushed into the grinding chamber from the front. The filter and the filter holder surface in the grinding chamber must be free from dirt and dust. Now the grinding chamber door can be closed and locked with the operating knob Pos. S 9. The screw mounting of the operating knob is screwed in by pressing and twisting it into the internal thread on the drive cover. The first turns are easy, if turning becomes difficult then the door may be closed, but the safety switch is not yet unlocked. Continue turning until stop (approximately two revolutions), only then is the safety switch unlocked and the mill ready to start.

Grinding operation with the cutting mill:

The suitable filter is pushed in and the door is locked before beginning the grinding procedure. Afterwards a collecting pan is fastened to the discharge cone Pos. S 1 (standard ground joint NS 29). Suitable containers are e.g. Erlenmeyer flasks Pos. Z 3, round-bottomed flasks, test tubes made of either glass or plastic. The containers are to be secured with a fixing clip Pos. Z 2 (see fig. 2).

For heat sensitive materials or with ground materials that have a very high increase in volume during cut up, it is advisable to use a two-necked flask with filter attached to its second (e.g. cotton wool, etc.). The ground material is cooled from the air flow and no excess pressure is created. For very critical ground materials (very fine powder after being cut up) the ground material can also be sucked off by means of a commercial vacuum cleaner with a suction tube adapted at the discharge cone. The powder is then in the vacuum cleaner bag where it can be removed for further processing steps.

The speed knob Pos. A 12 is turned to the left stop (3000 rpm). The mill is started by shifting the rocker switch Pos. A 13 to Position 1. First the green „Power“ LED is lit and then „Overload“ LED lights up briefly before the motor smoothly starts.

The speed can be set by turning the speed knob to the desired speed. A speed up to 5000 rpm is recommended for cutting-grinding while higher speeds should only be selected in exceptional cases (e.g. greater energy requirements such as for dry pieces of wood). The optimal speed must be determined by way of experiment. The higher the speed, the greater the cutting-grinding power. However, more energy is introduced into the ground material and thus the thermal load is increased.

The plunger Pos. 5 6 has to be removed from the guide channel before the ground material is fed in through the funnel Pos. 5 5. Make certain that the tubular side valve Pos. 5 8 is closed. Bulky materials can be simply filled into the funnel.

Metering is carried out by lifting the tubular slide valve.

Bulky material should be fed into the guide channel with the tubular side valve closed and the plunger must be fitted to in front of the tubular side valve in order to prevent chunks of ground material from spurting out.

The plunger may only be pulled out when the motor is running if the tubular side valve is closed.

Caution!

Note!

Caution: Wear protective personal equipment (protective glass, etc.)

You can hear if ground material gets into the grinding chamber and if too much material has been metered then the rotor speed decreases. This drop in speed is clearly indicated by the lower grinding noises and the orange „Overload“ LED lights up. Metering of the ground material should only be resumed when the speed has returned to the adjusted value or if the „Overload“ LED is no longer lit.

If this done then the mill becomes overloaded and switches off after some time because of overheating. This is then indicated by the red „Error“ LED. In addition, a beep sounds to indicate That the mill should be switched off using the rocker switch Pos. A13. The mill must then cool off for a period in order to be able to start it again (see chapter "Eliminating malfunctions").

When it comes to cleaning the grinding chamber experience has shown that it makes sense to

insert the plunger into the guide channel as far as possible and then to briefly switch the unit to maximum speed after the actual grinding operation has been completed and only then to turn the speed knob back to the left stop. The mill is then shut off with the switch.

Bulky and fibrous ground material, e.g. straw, hay, etc. should not be fed into the funnel lengthwise or in the form of larger positions. The plunger should be used to introduce them into the grinding chamber. Note that the plunger is coded and thus only fits into the guide channel in one direction (edge bevelling on the plunger and in the guide channel have to match each other). The „V“ on the plunger must be visible from the front. It has proved to be advantageous in some cases if the plunger is turned around and the bulky material is pushed from the funnel into the guide channel with the cylindrical end of the plunger.

The better a sample is prepared (pre-crushed to the optimal feed size), the easier the metering and the faster the grinding operation. In addition, failures such as rotor blocking can be virtually eliminated.

The filter is removed after opening the grinding chamber door and any remaining ground material is wiped into the collecting pan.

Note: After longer grinding cycles the grinding chamber and the filter may have an increased temperature. In this case wait until the parts have cooled off.

It may help to crush the ground material in several steps, i.e. with large perforations is used for pre-crushing and finer ones in each of the following steps. (see fig.4).

Advantages:

- Fewer failures, e.g. due to clogged filters or overloading the motor.
- Less heat-up of ground material since the retention time of the ground material in the grinding chamber is substantially reduced.
- Closer grain size spectrum for the final material.
- Normally less total grinding with better results as despite several runs comminution is performed more quickly and more effectively.

Disassembling the cutting-grinding head:

Caution: The blades on the rotor and in the grinding chamber have sharp edges. Observe safety regulations.

Before disassembling the head of the mill, the mill itself must first cool down as the grinding chamber and the beater rotor cannot be detached otherwise.

The power supply plug must be disconnected before disassembly.

After door is opened the countersunk screw Pos. S 21 is loosened by means of Allen wrench. After opening the door, the countersunk screw, P 15, is loosened using an Allen key, while the rotor is held

by the face spanner that is included in delivery. The spanner is positioned so that the pins of the spanner are inserted into the bore holes on the face of the rotor. The disc must be turned so that the bore hole in the rotor is accessible for the spanner pins.

If the hexagon socket of the screw is dirty then it must be cleaned by means of a brush. After the cheese head screw and the disc, S 22, have been removed the blade rotor, S 4, is removed. The adjustment spring must be in the area of the opening (in the back wall of the mill chamber). Now the three screw mountings Pos. 12 are loosened by means of a fork wrench and the grinding chamber is rotated in a clockwise direction and taken off.

Eliminating malfunctions

Mill does not start:

Possible causes:

- Power supply plug is not correctly plugged into the equipment and/or the socket („Power“ LED does not light up if the rocker switch Pos. A 13 is moved to Position 1).
- Grinding Chamber door not properly locked, („Power“ LED, „Error“ LED are lit and beep is heard when the rocker switch Pos. A 13 is moved to Position 1).

Remedy:

First shift the rocker switch Pos. A 13 to 0. Check voltage supply and the power supply plug is necessary. Or if „Error“ LED is lit and beeping can be heard close the grinding chamber door and tighten the turning handle Pos. P 10/ S 9 until stop. The rocker switch should only then be shifted to Position 1 again.

Engine hums after the rocker switch Pos. A 13 is shifted but the rotor does not begin to run.

Possible causes

- Ground material chunks is in the grinding chamber before grinding begins and blocks the rotor.

Remedy:

Switch the mill off at the rockers switch Pos. A 13 and disconnect the power supply plug. Open

the door and remove the ground material chunks from the grinding chamber. After the door is closed and the mill is connected to the mains again, you can switch the mill on again using the rocker switch Pos. A 13.

Engine switches off during operation:

Possible causes:

- Overloading of the engine because of too large feeding amounts which leads to heating up of the engine. Thus the motor thermal protective switch is triggered and switches the mill off. The overheating of the motor is indicated by the red „Error“ LED becoming lit and a beep is heard as well.
- Rotor blocking as a result of feed chunks that are too large or too tough.
- Rotor blocking due to improper speed adjustment (in general, the speed is too low for the given application).

Remedy:

Switch the mill off at the rocker switch Pos. A13 and disconnect the power supply plug. Open the door and remove any blocked ground material or ground material that has been cut up.

If the disturbance occurred because of overloading, let the mill cool off before you start it again.

Maintenance and cleaning

The grinding chamber interior (funnels, dosing equipment, rotor, grinding surface, filter and outlet nozzle) with a brush (no wire brush) should be cleaned after the grinding procedure. You can also use rags or cleaning brushes. The grinding surface does not need to be disassembled for cleaning between grinding operations with the same feed material. If a new grinding procedure with other feedings is to be performed whereby no contamination may be produced by the previously ground material, then more thorough cleaning is necessary which means dismantling the grinding chamber and the rotors. See the chapter „Commissioning“ for assembly and disassembly.

The grinding chamber parts can be washed with water and, if required, with a conventional dishwashing agent.

Be certain that the parts are dry before being re-assembled.

Caution!

Caution: Refer to the chapter on „Safety Instructions“ when cleaning the cutting-grinding surface and the blade rotor as the screwed on blades are very sharp and may result in injuries if they are improperly handled.

Clogged filters can be cleaned with brushes or compressed air. If ground material parts have settled into the filter holes in such a way that cleaning alternatives are not sufficient, then you will have to clean the clogged holes with suitable needles or clean the filter in an ultrasonic bath.

If beaters or blades and the grinding surface become worn in the course of the time, particularly with constant grinding of very hard substances, so that quick operation with the appropriate quality is no longer possible, then the worn out parts should be replaced.

Blades may not be reground as the blade clearance will become too large. For reasons of safety, only new, original IKA® bolts may be used for fastening the knives.

The appliance is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see **www.ika.com**, spare parts diagram and spare parts list.

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Only clean **IKA®** appliances using these **IKA®** approved cleaning agents. To remove use:

Dyes	Isopropanol
Building materials	Water containing detergent, Isopropanol
Cosmetics	Water containing detergent, Isopropanol
Food	Water containing detergent,
Fuels	Water containing detergent,

- Wear protective gloves when cleaning the devices.
- Do not place electrical appliances into the cleaning agents for cleaning purposes.
- Do not allow moisture to get into the appliance when cleaning.
- Please consult **IKA®** before using any cleaning or decontamination methods, other than those recommended here.

Accessories

Sieve insert with hole sizes:

Ø 0.25 mm
Ø 0.50 mm
Ø 1.00 mm
Ø 2.00 mm
Ø 3.00 mm

MF 10.1
MF 10.2

Cutting mill head
Impact grinding head

Technical data

Design voltage:	VAC	230±10%
or	VAC	115±10%
Design frequency	Hz	50/60
Power input	W	1000
Rotor speed	rpm	3000 - 6500
Peripheral speed impact grinding	m/s	35
Cutting speed cutting grinding	m/s	24
Cyclic duration	min	120/30
Intermittent operation		
Perm. ambient temperature	°C	+ 5 bis + 40
Perm. relative humidity	%	80
Protection class acc. DIN EN 60529		IP 22
Protection at overloaded		Current limitation
Fuses on mains plate		
230 V	A	6,3 Id.-No. 14 080 00
115 V	A	12 Id.-No. 26 493 00
Safety		Temperature switch on motor, door limit switch
Protection class		I
Overvoltage category		II
Contamination level		2
Noise level	dbA	max. 85
Operation at a terrestrial altitude	m	max. 2000
Dimensions (WxDxH)	mm	310x570x370
Weight	kg	10

Subject to technical changes!

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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GmbH & Co.KG