Powerful

Stirring

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

Operating Manual

Overhead Stirrers Hei-TORQUE

Operating manual must be read before initial start-up.
Please follow the safety instructions provided.
Please keep for future reference.

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44

About this Document	4
Version and variants	4
About this manual	4
Terms and signs used.....	5
Safety Instructions	6
General safety instructions.....	6
EC Declaration of Conformity.....	6
Intended use.....	6
Installation / Electrical safety	7
Personnel qualification.....	7
Operating company's obligations.....	7
Installation site	7
Changes to the unit	8
Safety for the personnel	8
Safety during use	8
Disposal	8
Device Description	9
Device overview	9
Hei-TORQUE Value overall view.....	9
Hei-TORQUE Value control panel	9
Hei-TORQUE Precision overall view	10
Hei-TORQUE Precision control panel	10
Start Up	11
Connect / disconnect power cord	11
Switch unit on and off	11
Operation	12
Control knob	12
Control keys	13
Hei-TORQUE Precision:	
Additional Functions.....	14
Symbols used.....	14
Menu navigation	15
Settings	16
Edit Profile.....	17
Favorites	19
Operation	21
Set Gear stage	21
Set speed.....	22

Rotation	23
Switch rotation on	23
Switch rotation off	23
Set rotation time	24
Start rotation via Favorites	24
Set rotation start time	26
Display progression graph	26
Calibrate Torque	27
Attach impellers	28
Hei-TORQUE Value:	28
Hei-TORQUE Precision:	30
Cleaning and Maintenance	31
Cleaning	31
Maintenance	31
Troubleshooting	32
Malfunctions and clear up	32
Assembly	33
Quick chuck	33
Support rod	34
Shaft guard	35
Electrical connections	36
Interface	37
Disassembly and Storage	37
Dismantling, transportation and storage	37
Dismantling	37
Transportation and storage	37
Accessories and Spare Parts	38
Scope of delivery	38
Accessories	38
Attachments	39
Technical data	39
Performance range	40
Service	41
Contact / Technical Service	41
Warranty	42
Confirmation of condition of unit	43

> Version and variants

Version

This manual describes the function, operation and maintenance of the stirrer Hei-TORQUE.

Version	Alteration date
1.1	04/2015

Variants

The units are available in different configurations. Certain features and functions apply to specific configurations only. The varieties are documented in this manual.

> About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

[http://www.heidolph-instruments.com/Support/Operation Manuals/Hei-Torque](http://www.heidolph-instruments.com/Support/Operation%20Manuals/Hei-Torque)

In addition please observe the regional regulations.

➤ Terms and signs used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the appliance.

The following terms and type of signs (basic drawing elements) are used:

Symbol	Additional term / Description
Warning signs 	<p>The yellow triangle indicates hazardous situations. It is used in combination with the following terms:</p> <p>DANGER: Indicates a hazardous situation which, if not avoided, will result in irreversible injury or death.</p> <p>WARNING: Indicates a hazardous situation which, if not avoided, could result in serious injury or material damage.</p> <p>CAUTION: Indicates a hazardous situation which, if not avoided, may result in a minor injury or material damage.</p>
Prohibitory signs 	<p>Not to do:</p> <p>The red circle indicates a situation that should be avoided under all circumstances to avoid injury or damage.</p>
Mandatory signs 	<p>To do:</p> <p>The blue circle indicates important information. Please obey to avoid property damage.</p>

Other signs used

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a. b. c.	List of variants

➤ General safety instructions

EC Declaration of Conformity

This device complies with the following EC-Directives*:

- 2006/42/EC Machinery Directive
- 2014/35/EU Low-Voltage Directive
- 2014/30/EU Electromagnetic Compatibility Directive

*Also see attachment "EC Declaration of Conformity".

The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure the operating manual is available at all times.

The device may only be used under the following circumstances:

- Only operate the device, if it is in full working order.
- Ensure all operators of the device possess the necessary safety and risk awareness.
- Operate the device according to the instructions stipulated in this manual only.
- If there is something you do not understand, or certain information is missing, ask your manager or contact the manufacturer.
- Do not do anything on the device without authorization.
- Only use the device according to its intended use.

Intended use

The device is intended for use by trained and authorized personnel only.

The device is suitable for the following use:

- stirring
- mixing
- blending of solvents and/or solids

The operation for the device is suitable in one of the following locations:

- in any research application

Improper use

Any use which deviates from the device's intended use is considered to be improper. The manufacturer does not accept liability for any damages resulting from non-permitted uses. The risk is carried by the operator alone.

Installation / Electrical safety

- The device may only be connected when the mains voltage corresponds to the information on the type plate of the unit.
- The mains connection must be accessible at all times.
- Repairs may only be performed by a qualified electrician.
- Never operate the unit with a damaged power cord.
- Always turn the unit OFF and disconnect mains power before performing any maintenance or service.

Personnel qualification

- The device may only be operated by qualified persons.
- The device may only be operated by individuals who have been instructed in its proper use by qualified persons.
- The device may only be operated and maintained by persons who are of legal age.
- Other personnel may only work with the unit under continuous supervision of an experienced qualified operator.
- This manual must be read and understood by all persons working with the device.
- The personnel must have received special safety instructions in order to guarantee responsible and safe work procedures.

➤ Operating company's obligations

Installation site

- The device must be positioned in a suitable location.
- The device must be installed sufficiently stable.
- All screw connections must be securely tightened.
- Installation and operation of the device is only permitted in facilities which are fitted with the appropriate laboratory equipment (e.g. with air extraction units).
- We recommend operation in fume hood enclosures.
- The device may only be operated in enclosed spaces and under the following environmental conditions:

Ambient temperature	<ul style="list-style-type: none">▪ 5 - 31 °C at 80 % relative humidity, no condensation▪ 32 - 40 °C decrease linearly till max. 50% relative humidity▪ up to 2000 m (6500 feet) height above sea level
---------------------	---

Depending on the medium used:

- ➔ Only operate the device in conjunction with an extractor hood (at least 10-fold air change, with error monitoring), see DIN EN 14175 and DIN 12924.

Changes to the unit

- No unauthorized changes may be made to the unit.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the appliance may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the unit other than described in this manual.

Safety for the personnel

- Ensure that only qualified personnel operate the device.
- Observe the following regulations:
 - Laboratory guidelines
 - Accident prevention regulation
 - Ordinance on Hazardous Substances
 - Other generally accepted rules of safety engineering and occupational health
 - Local regulations

➤ Safety during use

- Wear the appropriate protective clothing when working on the device (protective glasses and if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic waves) for the frequency range (3×10^{11} Hz to 3×10^{15} Hz).
- Avoid putting pressure on the display when you do not operate the device.
- Eliminate errors immediately.
- Do not use abrasive material to clean the glass surfaces. Only wipe with damp cloths.
- Always switch the mains off after use.

➤ Disposal

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device according to the valid national legal regulations.
- Dispose of the packaging material in accordance to the valid national legal regulations.


➤ Device overview

Hei-TORQUE Value overall view


Hei-TORQUE Value control panel


Device Description

Hei-TORQUE Precision overall view

Hei-TORQUE Precision control panel


➤ Connect / disconnect power cord

Connect power cord

The plug comes with a yellow clip on the bottom side to ensure a firm attachment.

Power cord

Clip

Plug

- Connect the plug to the socket on the back side of the unit with the yellow clip pointing downwards.
- ✓ The plug is securely attached once the yellow clip clicks in place.
- Connect the power cord on the other side to the electrical ground socket.

Socket

Disconnect power cord

- Press the yellow clip on the downside of the plug and hold it in this position.
- Pull the plug simultaneously out of socket of the unit.

➤ Switch unit on and off

Switch unit on

- Press the on/off button.
- ✓ The button locks in.
- ✓ The display is illuminated.
- ✓ The screen displays a short animation on how the rotation is switched on.
- ✓ The according touch panel flashes for 2 seconds.

Switch unit off

- Press the on/off button.
- ✓ The display powers off.

On/Off switch

➤ Operation

Control knob

The knob is attached to the panel magnetically. Settings can be adjusted without the knob by rotating your finger on the touch screen.


As soon as the knob is removed a diagram is displayed for 3 seconds showing how to change parameters without the knob.

- Adjust your settings by rotating your finger clockwise to increase or counter clockwise to decrease the settings.

Hei-TORQUE Value: speed adjustment

- Turn control knob clockwise.
- ✓ The speed increases.
- Turn control knob counter-clockwise.
- ✓ The speed decreases.

Hei-TORQUE Precision: speed adjustment

- to address and change settings: turn control knob
- to select and confirm settings: push control knob
- Turn control knob clockwise.
- ✓ The indicator jumps to "Timer" or "Speed".
- Push control knob to confirm selection.
- ✓ Set figures change from white to orange for approximately 2 seconds.
- While figures are orange turn control knob clockwise to increase or counter clockwise to decrease the setting.
- Push control knob to confirm setting.

Start Up

Settings can only be changed while figures are orange.

Control keys

Hei-TORQUE Value 400

Speed gear is set via the touch-sensitive control key.

- Touch control key for gear stage.
- ✓ The operating panel will display the actual gear setting.

Control key gear stage


Hei-TORQUE Precision

The unit is equipped with several touch-sensitive control keys.

The functions mentioned above the control keys apply when touching the control key.

- ✓ In the main menu adjacent three control keys are touch-sensitive.
- Touch one of the control keys.
- ✓ The above mentioned function applies.
- ✓ In the submenu adjacent two control keys are touch-sensitive.
- ✓ The control key "back"  leads you back to the previous step.
- ✓ The control key "home"  leads you back to the main menu.

➤ Hei-TORQUE Precision: Additional Functions

Symbols used


Symbol	Naming	Function
	Home	Back to main menu
	Back	Back to previous menu
	Save	Save setting or process
	Ramp jump	Rotation speed jumps immediately to predefined setting
	Ramp slope	Rotation speed increases or decreases gradually to predefined setting
	Correction	Deletes last input
	Delete	Deletes template
	Continuous operation	Activates repetition of process cycle
	Activate	<ul style="list-style-type: none"> ▪ Activates single operation of profile ▪ Activates input
	Deactivate	<ul style="list-style-type: none"> ▪ Deactivates profile ▪ Deactivates input
	Profile cycle	Profile cycle active: Adjustment of time or speed not possible
	Continuous operation	Continuous operation active: Adjustment of time or speed not possible
	Start time	Start time active

Menu navigation


Switch to submenu

→ Touch control key „Menu“.


Control key Menu

- ✓ The submenu opens.
- Turn control knob clockwise or counter clockwise.
- ✓ The position indicator jumps to a menu item.
- ✓ The color of addressed menu item changes from white to orange.

- Push control knob.
- The selected menu will be displayed.

Quit sub menu

- Touch control key “Back”
- ✓ You will return to the previous menu.
- Touch control key “Home”
- ✓ You return to the main menu.

Settings

You can set the time and language through the menu item "settings"


Set language

- Open the submenu "settings".
- Select menu item "language" and choose your desired language.
- Push control knob to confirm selection.
- Leave the menu "settings" with the control key "back" .

Set time

- Select menu item "time".
- Push control knob to adjust time.
- Turn control knob to change time.
- Push control knob to complete setting.
- Push control knob again to confirm and save selection with "OK".
- Leave the menu "settings" with the control key "back" .

Start Up


Save Profile

All profile steps are done:

- Turn the control knob left until the saving symbol is assigned in orange color.
- Push the control knob to confirm selection.

- ✓ An alphanumeric block appears to enter the name of your Favorite.
- Choose a letter or number by turning the control knob.
- Push the control knob to confirm selection.
- Delete accidentally inputs through the deletion symbol .
- Save the selection with "OK".
- ✓ The confirmation "data saved" appears shortly.
- Leave the menu with to enter a further Favorite.

- ✓ The last saved profile appears as a template.
- Adjust settings and save profile with a new name.

Choose your name accurately since it's possible to assign a name twice .

Change speed gear for the next profile:

- ✓ The last saved profile appears as a template.
- Choose and confirm the delete-symbol .
- ✓ The last saved profile step is deleted from the template and you receive a free menu.
- Leave the menu with the control key "home" .
- Choose a different speed gear in the main menu.
- Return to the menu "Edit Profile".
- Enter a new profile with changed speed gear.

- One profile can only be saved with either speed gear I or speed gear II.
- The speed gear for the profile must be chosen in the main menu.
- Is the Template for a new profile in speed gear I (or II), the new process has to be entered in speed gear I (or II).


Favorites

Display Favorites

- Open the menu item "Favorites" in the submenu.
- ✓ Saved favorites from menu item "Profile edit." are shown.
- Choose a favorite and confirm the selection.

- You receive adjacent selection.
- Choose "Show".
- The saved profile is shown.

The symbols for "single operation" and "continuous operation" are displayed on the bottom part of the display. You can select one of the symbols however you will return to the previous profile.

Choose a Favorite for a Template

If you want to use a particular Favorite for a template, don't go to the menu "Edit Profile" but to the menu "Favorites". Save your Favorite with a new name under the menu "Edit".

Start Up

Edit Favorites

Choose from following options:

- Open an existing profile as a template for a new profile
 - Delete profile templates to edit a new profile
 - Start a profile
- Open the Favorite you want to edit
 - Choose "edit".
 - ✓ The saved Profile appears in your Favorites.

Edit


- Choose the row you want to edit with the cursor.
- Push the control knob to release the row for editing.
- Change the desired data.
- Save the changes with  under a new name.

Delete template

- Klick the delete symbol .
- The template except the edit line is deleted.

Start

- Klick the symbol "single operation"  or "continuous operation" .
- ✓ The symbols appear on the upper part of the display.

An already saved Favorite can't be changed supplementary.

If you want to change an existing Profile, open the Profile and save it under a new name.

➤ Set Gear stage

Only available for Hei-TORQUE Value 400 und Hei-TORQUE Precision 400

To ensure constant speed under changing loads you can choose between two gear settings.

Speed range:

- Gear stage I : 10 rpm – 400 rpm
- Gear stage II : 200 rpm – 2000 rpm


Gear transmission can only be changed, if the unit is not in operation.

- ✓ Rotation is switched off.
- ➔ Touch the gear stage on the operating panel.
- ✓ The operating panel will display the actual gear setting.
- ✓ Speed settings will automatically adjust to the lowest rpm.

Hei-TORQUE
Value 400

Hei-TORQUE
Precision 400

Control key gear stage

Hei-TORQUE Precision 400

The speed gear does not switch automatically when choosing a Favorite.

Always set the suitable speed gear in the main menu before you choose a Favorite to start a Profile.

➤ Set speed


The overhead stirrer operates clockwise. The speed range varies from

- 10 rpm - 2,000 rpm

(For more details see chapter „Technical data“.)

Hei-TORQUE Value: set speed

- ➔ Turn control knob clockwise.
- ✓ The speed increases.
- ✓ The speed setting is documented on the display under „set rotation speed“.
- ➔ Turn control knob counter-clockwise.
- ✓ The speed decreases.

Increase rotation

Decrease rotation

Hei-TORQUE Precision: set speed

- ✓ You are currently in the main menu.
- ➔ Turn the control knob clockwise to the right until the cursor is located left to "speed".
- ➔ Push the control knob.
- ✓ The color of set-rpm changes 2 sec. long from white to orange.
- ✓ Turn the control knob clockwise to the right or left to increase or decrease speed.

Hei-TORQUE Value: adjust speed

At any time during operation turn the control knob clockwise or counter-clockwise to adjust speed.

Hei-TORQUE Precision: adjust speed

- ✓ You are currently in the main menu and have set and start speed.
- ➔ Choose "speed" .
- ➔ Increase or decrease speed.

The speed cannot be changed if the Profile was already activated through Favorites

Rotation

Warning: **Risk of injury!**

Rotating items may partly injure the operator if in contact with hair, clothing, jewelry and other items.

- Never touch any rotating parts.
- Don't wear clothing belts, ties which may get gripped by rotating parts.
- Don't wear long necklace or bracelets.
- To protect long hair please wear a cap.

Warning: **Risk of poisoning!**

Depending on your application stirring may result in a chemical reaction which sets free dangerous gas.

Operate the unit eventually behind a closed and ventilated fume hood.

Depending on viscosity, speed and volume media may splash after switching on the unit.

Initially choose a low speed range and increase rpm step by step.


Switch rotation on

- Touch the slider with one finger and move from left  to right .
- ✓ The unit starts rotation.
- ✓ Rotation speed ramps up in a short period until it has reached the pre-programmed rpm-settings.

Switch rotation off

- Touch the slider with one finger and move from right  to left .
- ✓ Rotation stops.

Hei-TORQUE Precision

You can also stop rotation via slide panel, if you have activated a Profile through Favorites .

Set rotation time

- The adjustable time suffice from 1 minute to 24 hours.
- The sequence is displayed to the split second.

- Turn the control knob clockwise to the right.
- ✓ The positioning indicator jumps left next to "Timer".
- Push the control knob.
- ✓ The color of set-rpm changes 2 sec. long from white to orange.
- Turn the control knob clockwise to the right or left to increase or decrease time.
- Start the rotation via slide panel.
- ✓ The scale below set-time indicates the remaining time.
- ✓ In spite of the set-time, the rotation can be stopped via slide panel at any time.

Start rotation via Favorites

There are two ways to start a saved profile:

- Single operation
 - 
 - The profile operates one-time on basis of the pre-programmed time.
 - The rotation stops automatically as soon as the operation is finished.
 - The profile symbol appears on the upper part of the display.
- Continuous run
 - 
 - The profile operates until the end.
 - After the cycle is finished, the operation starts again.
 - The operating continuous until the profile is stopped manually.

Operation

Single operation

- Open the Favorite you want to activate.
- Choose "Edit".
- ✓ The deposited profile is displayed.
- Choose "single operation"
- ✓ The symbol "profile active" appears on the upper part of the display.
- Return to the main menu.
- Start the profile via slide panel.
- ✓ The profile operates and ends after one-time.
- Stop the profile through the menu "Favorites", "Edit" by deactivation
- ✓ The symbol for "profile active" disappears.

Continuous operation

- Open the Favorite you want to activate.
- Choose "edit".
- ✓ The deposited profile is displayed.
- Choose "continuous operation"
- ✓ The symbol "profile active" and the symbol for "continuous operation" appear on the upper part of the display.
- Return to the main menu.
- Start the profile via slide panel.
- ✓ The profile starts operation and continues.
- ✓ Stop the continuous operation via slide panel.
- Stop the profile in the menu "Favorites", "Edit" by deactivation

Set rotation start time


Set a start time:

- for one of your Favorites.
- for a setting which you have directly conducted in the main menu. In due consideration of the pre-programmed time.

The start time is interpreted as the actual time. To assure an accurate start it's important to set the clock correctly to the time of day.

- Activate the profile for one of your Favorites or conduct the settings in the main menu.
- Open the menu item "start time" in the submenu.
- Push the control knob to edit the time.
- Set the time for profile start.
- Push control knob to confirm set time.
- Activate the set time with "activate" .
- ✓ The rotation starts to the pre-programmed time.

Display progression graph

It's possible to display the profile, on the basis of a diagram, while operating a profile. Monitor the progression especially for samples with high viscosity and for edit profiles.

The diagram is only displayed when rotation is active.

Operation



Hei-TORQUE Precision 100 / 200

- You have currently started stirring.
- Choose the menu "Graph".

Hei-TORQUE Precision 400

- You have currently started stirring.
- Open the submenu.
- Choose the menu item "Graph".

- ✓ A diagram of the current profile appears.
- ✓ The diagram shows the last 3,5 minutes of the profile.

Calibrate Torque

It's reasonable to assign the Torque individually for samples with a different start and end viscosity.

- The change from higher to lower viscosity of samples is indicated by a negative torque.
- The change from lower to higher viscosity of samples is indicated by a positive torque.
- ✓ You have currently started stirring.
- Choose the control key "Cal" in the main menu.
- ✓ The notification "Torque calibrated" appears.
- ✓ The display shows positive or negative torque, depending on viscosity of your sample.

Reset calibration

You can reset calibration any time with the control key "cal" during operation as well as standstill.

➤ Attach impellers

Warning:
Risk of injury!

If the unit is not powered off completely, rotation may start unexpectedly while opening or closing the quick chuck. This may result in serious injuries.

Always switch the unit off when opening and closing the quick chuck.

Hei-TORQUE Value:

Open quick chuck

The device is switched off.

1. If an impeller is already fixed:
 - Hold the orange ring on the quick chuck and turn the black ring clockwise until the quick chuck opens.
 - Hold the impeller with one hand and turn the black ring further clockwise until the impeller is completely disconnected from the unit.
2. If no impeller is fixed:
 - Turn the black ring clockwise until the rod can be moved through the quick chuck.

Quick chuck

Impeller

Operation

Close quick chuck

**Warning:
Risk of injury!**

If the quick chuck is not closed properly the impeller may disconnect during rotation and result in serious injuries.

Impeller is securely attached once the quick chuck has locked in.

- Move the impeller upwards until it has reached an operational position.
- Hold the impeller with one hand and turn the black ring counter-clockwise until the quick chuck locks in.

Quick chuck

Impeller

Impeller is securely attached once the quick chuck has locked in.

Do not close the black ring further. Otherwise the quick chuck will be hard to open again.

If accidentally the quick chuck was closed to tight to open, demount the complete quick chuck. The quick chuck will be handled easier (refer to chapter "Assembly - Quick chuck" and "Disassembly - Dismantling").

Hei-TORQUE Precision:

Open Quick chuck

- The unit is switched off.
- Turn the locking ring clockwise 90° to the right.


✓ The locking ring is open.

- Pull down locking ring so that the orange ring is covered.
- Turn quick chuck counter-clockwise to the left beyond the resistance to open.
- Remove impeller.

Install impeller

- Insert impeller to the quick chuck.
- Close quick chuck carefully until it encounters resistance of the impeller.
- Adjust impeller axial.
- Close quick chuck beyond the first resistance until it clicks into place.
- Move up locking ring as far as it will go.
- Turn locking ring counter-clockwise to the left as far as it will go.

✓ The impeller is installed.

› Cleaning

The quick chuck and all surfaces of the unit allow for cleaning with a damp cloth and if required a mild soap lotion.

Attention:
Damage to the surfaces

All surfaces can be damaged by improper cleaning.

In any case do not use:

- Chlorine bleach or any cleaning agent including chlorine
- Ammonia
- Abrasive cleaning agents such as cleaning rags, scrubbing agents or any other agents which include metal components

› Maintenance

The unit is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

Troubleshooting

➤ Malfunctions and clear up

Hei-TORQUE Value

Malfunctions	Possible reason	Troubleshooting
Display cannot be switched on	No power	<ul style="list-style-type: none"> Check size and shape of plug and compatibility to your electrical socket. Check circuit breakers
Rotation stops and torque display blinks	Overload protection has switched the unit off	<ul style="list-style-type: none"> Switch off the device and allow it to cool down. Reduce rotation speed. Use a different impeller.

Hei-TORQUE Precision

Malfunctions	Possible reason	Troubleshooting
Display cannot be switched on	No power	<ul style="list-style-type: none"> Check size and shape of plug and compatibility to your electrical socket. Check circuit breakers
No Rotation	Symbol  at the upper part of the display: profile active	Stop profile through the menu "Favorites" and "edit" or "show".
Error: „Quick chuck not fixed“	Quick chuck closed incorrectly	Close Quick chuck (see chapter "Install impeller", section "Hei-Torque Precision").
Favorite won't reach selected set rpm	Wrong gear stage active	<ul style="list-style-type: none"> Stop process. Change gear stage. Start process.
Rotation stops and torque display blinks	Overload protection has switched the unit off	<ul style="list-style-type: none"> Switch off the device and allow it to cool down. Reduce rotation speed. Use a different impeller.

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph Instruments distributor immediately.


Quick chuck

With delivery the following pieces are provided:

- a quick chuck to hold the impeller
- an Allen-key to install the quick chuck

Mount Quick chuck

- Push Quick chuck onto the retaining bolt by sliding tongue and groove into each other.
- If it does not move smoothly, turn the allen screw counter-clockwise approximately 10° out of the quick chuck and retry.
- Tighten Allen-screw through the drilling to fix the quick chuck.

Hei-TORQUE Precision

To positioning the hexagon socket screw key, turn the supporting bolt with the orange ring counter-clockwise to the left until the hexagon socket screw key is visible below the drilling.

The hexagon screw is very small and can be easily lost. If you want to demount the quick chuck, loosen the hexagon screw only approximately 60° to prevent it from accidentally falling out.

➤ Support rod

On the back side of the unit you will find a slot to install the support rod.


To install the support rod use the supplied TORX L-key.

The back support rod allows you to attach the stirrer to a stand in conjunction with a 2-way clamp or to a wall grid.

→ Install the support rod in the slot on the back side of the unit.

→ Ensure that the threaded hole in the support rod aligns with the drill hole in the housing.

→ Secure the support rod with the threaded screw supplied.


➤ Shaft guard

To prevent injuries caused by rotating parts a shaft guard can be attached. (optional accessory, see section "Accessories")

Attach shaft guard

- Hold the adaptor to the quick chuck with the cut-out below the on/off switch.
- ✓ The two guide aids ensure correct fitting.
- Attach the adaptor with the two allen screws and key supplied to the bottom side of the housing and tighten them until it stops.
- Attach both knurled screws to the drill hole on the adaptor.
- Push the grooves on the shaft guard alongside the knurled screws until the screws reach the end of the groove.
- Tighten both knurled screws to hold the shaft guard firmly in place.

Remove shaft guard

- Untighten both knurled screws which hold the shaft guard to the adaptor.
- Turn the shaft guard slightly left and remove it from the stirrer.

Electrical connections

- This unit must **only** be connected **by a professional electrician**.
- **Repairs** of any kind are allowed by **qualified professionals only**. Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.


Warning: Electric shock!

This unit must only be connected to a grounded electrical socket.

The unit and voltage must match. The type label on the back side of the unit provides all voltage specifics (see picture below).

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

Color code for electrical connections:

Green/Yellow	PE: Protective conductor (Earth)
Blue	N: Neutral conductor
Brown	P: Phase
Black	P: Phase
Grey	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is approved by the local regulations.

➤ Interface

A mini USB-interface is located on the back of the unit.

You can connect the unit to a PC via the USB-interface.

USB interface

➤ Dismantling, transportation and storage

Dismantling

Warning:
Risk of injury!

While the unit is plugged in you are at harm if the rotation is accidentally switched to ON mode.

Prior to dismantling always disconnect the plug from the electrical socket.

- ✓ The unit is switched off.
- Disconnect plug from power socket and unit socket.
- Remove impeller and eventually dismantle shaft guard.
- Remove quick chuck.

Demount quick chuck

- Only Hei-TORQUE Precision: Position Allen-screw below Drilling.
- Hold on to quick chuck.
- Turn hexagon screw approximately 60° to 70° counter-clockwise to the left.
- ✓ The quick chuck can be pulled off.

Transportation and storage

- Store and transport the unit and its components only if they were emptied and cleaned in the original packing material.
- Alternatively store and transport the unit in a container which prevents damages during transportation.
- Seal the packing carefully against unauthorized or mistaken opening.
- Store the unit in a dry place.

Improper transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

➤ Scope of delivery

Component	Variant	Quantity	Product number
Hei-TORQUE	Value 100	1	501-61010-00
	Value 200	1	501-62010-00
	Value 400	1	501-64010-00
	Precision 100	1	501-61020-00
	Precision 200	1	501-62020-00
	Precision 400	1	501-64020-00
Quick chuck Hei-TORQUE Value		1	509-02300-00
Quick chuck Hei-TORQUE Precision		1	509 02200 00
Support rod		1	22-02-14-01-41
Rotation knob		1	11-300-007-19
Operating Manual English / German		1	01-005-005-55
Warranty registration / Confirmation of condition		1	01-006-002-78
EC Declaration of Conformity		1	01-001-025-01
Power cord		1	

➤ Accessories

Component	Quantity	Product number
Shaft guard	1	509-08100-00

Additional accessories are shown in the general catalog or on our website at www.heidolph-instruments.com

➤ Technical data

Hei-TORQUE

Standard supply voltage	1 x 230 V (50/60 Hz)
Operational voltage	85 - 260 V (50/60 Hz)
Protection class (DIN EN 60529)	IP54
Interference suppression	EN 50081
Drive	EC-Motor
Overheat protection	KTY-temperature control, self reset
Overload	automatic power shut off, message in display
Speed control	electronic
Noise emission (dB)	< 50
Shaft (Ø x L) (mm)	13 x 160
Ø Impeller min./max. (mm)	2.5 - 12.5
Permissible duty cycle	continuous
Permissible ambient conditions	<ul style="list-style-type: none"> ▪ 5 - 31 °C at 80 % relative humidity, no condensation ▪ 32 - 40 °C decrease linearly till max. 50% relative humidity ▪ up to 2000 m (6500 feet) height above sea level

Modell Hei-TORQUE

	Value	Precision
Interface	no	USB serial
Timer	no	yes
Control panel	TFT-Display, monochrom, 2.4"	TFT-Display, color, 3.2"
Torque	symbol	figure

Modell Hei-TORQUE

	Value/Precision 100	Value/Precision 200	Value/Precision 400
Supply power (W)	90	120	150
Stirring cap. H ₂ O max. (L)	50	50	100
Viscosity, max. (mPa s)	60,000	100,000	250,000
Torque, max. (Ncm)	100	200	400
Number of speed gears	1	1	2
Speed range (rpm)	10-2,000	10-2,000	Gear I 10-400 Gear II 200-2,000
Weight (g)	4,150	5,050	5,200


Modell Hei-TORQUE

		Value 100/200	Precision 100 /200	Value 400	Precision 400
Dimensions (mm)	B	86	86	93	93
	H	257	273	257	273
	T	241	241	241	241


Performance range


Hei-TORQUE Value 100
Hei TORQUE Precision 100

Hei-TORQUE Value 200
Hei TORQUE Precision 200

Hei-TORQUE Value 400
Hei TORQUE Precision 400

➤ Contact / Technical Service

Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local, authorized Heidolph Instruments distributor.

Warning: Danger of poisoning!

Contaminated units can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please

- advice exact substance
 - take adequate protective action towards our parts receiving and service personnel
 - mark the pack in accordance with Ordinance on Hazardous Substances
-

Included at the end of this manual you will find a „Confirmation of condition of unit“.

→ Kindly copy and fill in this form and submit it prior to shipping the unit for repair.

Contact details:

Heidolph Instruments Germany

Heidolph Instruments GmbH & Co. KG
Technical Service
Walpersdorfer Str. 12
D-91126 Schwabach / Germany

Tel.: +49 – 9122 - 9920-74
Fax: +49 – 9122 - 9920-84

E-Mail: service@heidolph.de
www.heidolph.de

Heidolph Instruments North America

Phone: 1-866-650-9604
E-mail: service@heidolph.com
www.heidolphNA.com

Heidolph Instruments United Kingdom

Phone: 01799 - 5133-20
E-mail: service@radleys.co.uk
www.heidolph-instruments.co.uk

All other countries

You will find contact details of your local Heidolph distributor at www.heidolph.com

➤ Warranty

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online (www.heidolph.com). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

Warranty claim?

- ➔ Please inform Heidolph Instruments or your local distributor should you wish to make a warranty claim.

➤ Confirmation of condition of unit

→ In the case of repair, copy and complete the Confirmation of condition of unit and send it to Heidolph Instruments.

1. Details about the unit

Product number _____

Serial number _____

Reason for repair _____

2. Has the device been cleaned, decontaminated/sterilized?

Yes _____ No _____

3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes _____ No _____

If not, which substances has the unit come into contact with?

4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date _____ **Signature** _____

Company stamp

Please note

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG