



**Exercise due care in working with solvents! Follow carefully all safety guidance in the current safety data sheet! Always operate the apparatus in a fume hood!**



**Glass can break and cause injury! In working with glass components, observe all appropriate safety precautions.**



**Danger of electric shock! Make sure that no liquids get into the cable connections or the inside of the equipment.**




**Caution: vessels get hot and can cause burns! Do not touch the glassware with bare hands during and immediately after an extraction.**

### List of Components

- Stand, complete with heating element, flask adapter and distance part
- behrotest® glassware kit, consisting of:
  - Round-bottomed flask with NS 29 std. taper neck, 250 ml (KEX 100) or 500 ml (KEX 250)
  - Soxhlet extractor EZ 100 / EZ 250
  - Reflux condenser RFK 100
- Extraction thimbles EX 100 HS or EX 250 HS (2 pcs)

### Setting up and Connecting

- Place the stand on a flat level surface.
- Remove the shipping protection from the heating unit.
- Verify that the power dial on the heating unit is set to the "0" off position.
- Insert the power plug into a mains socket.



212107

### Installing the Glassware

- 1 Insert the round-bottom flask into the heating unit.
  - 2 Insert the extractor in the mouth of the flask, threading it through the gaps in the mounting.
  - 3 Insert the reflux condenser on top of the extractor.
- In the bottom of the heating adapter you'll notice a cross-shaped profile. It serves as a distance piece to prevent the round-bottom flask from being jammed, and possibly crushed, in the adapter.

### Connecting the Tubing

The lower connection on the reflux condenser is the cooling water inlet and the upper connection is the water outlet.

Connect the tubing as follows:

- Slip a hose on the lower nozzle of the condenser and connect it to the laboratory water tap.
- Slip a hose on the upper nozzle of the condenser and lead it to a sink or drain.




### Turning the Heating Unit On

Switch the heating unit on by turning the thumb-wheel at the bottom of the unit.

To heat the solvent in the round-bottomed flask as rapidly as possible, it is advisable to initially set heating power to the uppermost level. Once the solvent starts boiling, the power level should be lowered.

### When the Extraction is Done...

... just place the reflux condenser in the boring in the stand so you have it handy for the next extraction.



212107

## Spare Parts

Art. No.	Designation	Description
<b>For 100-ml extraction (KEX 100):</b>		
B 0021 8499	RK 250	Round-bottomed flask, NS 29 neck, 250 ml
B 0021 7967	EZ 100	Extractor
B 0021 8214	RFK 100	Reflux condenser
B 0011 6487	EX 100 HS	Extraction thimbles, 33 x 94 mm, 25 pcs
<b>For 250-ml extraction (KEX 250):</b>		
B 0021 8500	RK 500	Round-bottomed flask, NS 29 neck, 500 ml
B 0021 7974	EZ 250	Extractor
B 0021 8214	RFK 100	Reflux condenser
B 0021 7975	EX 250 HS	Extraction thimbles, 33 x 205 mm, 25 pcs

## Technical Data

Nominal voltage	230 / 115 V~, 50 / 60 Hz
Current	at 230 V: 3 A; at 115 V: 6 A
Power consumption	450 W
Weight	approx. 7,5 kg
Dimensions (WxHxD, mm)	approx. 230 x 330 x 800

**behr Labor-Technik GmbH**  
 Spangerstraße 8  
 D-40599 Düsseldorf

Telephone: (+49 211) 7 48 47 17  
 Telefax: (+49 211) 7 48 47 48  
 eMail: info@behr-labor.com  
<http://www.behr-labor.com>

## Compact Apparatus KEX 100 / KEX 250

For 100-ml / 250-ml Soxhlet Extraction

