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CertoClav Vacuum Pro 12-22L

Laborautoklav

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Beispiel:

The screenshot shows the 'Settings' screen with a sidebar menu on the left containing: User Account, Network, Device, Audit Logs, Notifications, Sterilization, Calibration, and GLP. The main area lists sensor offsets: Steam Sensor (-1.0 °C), Heater Sensor (0.2 °C), Steam Generator Sensor (0.3 °C), Media Sensor (0.1 °C), and Pressure Sensor (-10 kPa / bar). The Pressure Sensor value is highlighted with a red box. Below the list is an 'Apply Offset' button. The bottom status bar displays 'STEAM: 89 °C', 'MEDIA: 86 °C', 'PRESS: 2 kPa' (highlighted with a red box), and 'support@certoclav.com'.

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This screenshot is identical to the previous one, but the Pressure Sensor offset is now set to -12 kPa / bar, which is highlighted with a red box. The 'Apply Offset' button remains visible. The status bar still shows 'PRESS: 2 kPa' in red.

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This screenshot shows the 'Settings' screen after the calibration. The Pressure Sensor offset is still -12 kPa / bar. The 'Apply Offset' button is now highlighted with a red box. The bottom status bar shows 'STEAM: 30 °C', 'MEDIA: 27 °C', 'PRESS: 0 kPa' (highlighted with a red box), and 'support@certoclav.com'.

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

Settings			
User Account	Time	User	Event ID
	Feb 03, 2021 02:45 PM	Admin	33
Network	Time	User	Event ID
	Feb 03, 2021 02:45 PM	Admin	35
Device	Time	User	Event ID
Audit Logs	Feb 03, 2021 02:44 PM	Admin	18
Notifications	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	18
Sterilization	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	18
Calibration	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	15
GLP	Time	User	Event ID
STEAM: 89 °C MEDIA: 86 °C PRESS: 1 kPa support@certoclav.com			

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Settings			
User Account	Time	User	Event ID
	Feb 03, 2021 02:45 PM	Admin	33
Network	Time	User	Event ID
	Feb 03, 2021 02:45 PM	Admin	35
Device	Time	User	Event ID
Audit Logs	Feb 03, 2021 02:44 PM	Admin	18
Notifications	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	18
Sterilization	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	18
Calibration	Time	User	Event ID
	Feb 03, 2021 02:44 PM	Admin	15
GLP	Time	User	Event ID
STEAM: 81 °C MEDIA: 78 °C PRESS: 1 kPa support@certoclav.com			

Time	User	Event ID
Feb 03, 2021 02:44 PM	Admin	18

Old Value	New Value
preferences_list_pressure_unit = bar	preferences_list_pressure_unit = kPa

OK

12.4.9 Audit Trail Event ID's

Event ID	Beschreibung
1	ACTION_OPEN_WINDOW
2	ACTION_USER_CREATED
3	ACTION_USER_CHANGED
4	ACTION_USER_DELETED
5	ACTION_PROGRAM_STARTED
6	ACTION_PROGRAM_INDICATOR_CHANGED
7	ACTION_PROGRAM_CANCELED
8	ACTION_PROGRAM_FINISHED
9	ACTION_PROGRAM_CREATED
10	ACTION_PROGRAM_CREATED_FAILED
11	ACTION_PROGRAM_DELETED
12	ACTION_PROGRAM_DELETION_FAILED
13	ACTION_PROGRAM_SAVED
14	ACTION_PROGRAM_SAVE_FAILED
15	ACTION_PROGRAM_FAILED
16	ACTION_SUCCESS_LOGIN
17	ACTION_FAILED_LOGIN
18	ACTION_LOGOUT
19	ACTION_ACTIVATED
20	ACTION_PREFERENCE_SETTINGS_CHANGED
21	ACTION_ERROR

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

17.3 Anhang 3: EMV

Electromagnetische Emissionen

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Steam sterilizer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Vacuum	The Steam sterilizer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	2 kV for power supply lines 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	1 kV line(s) to line(s) 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Steam sterilizer requires continued operation during power

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
IEC 61000-4-11	(60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	(60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	mains interruptions, it is recommended that the Steam sterilizer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 150 kHz to 80 MHz 3 Vrms 3 V/m 80 MHz to 2,5 GHz 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Steam sterilizer including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2,3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Steam sterilizer is used exceeds the applicable RF compliance level above, the Steam sterilizer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Steam sterilizer. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distances between

Portable and mobile RF communications equipment and the Steam sterilizer

The Steam sterilizer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Steam sterilizer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Steam sterilizer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
		m	
W	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = 1,2$	$d = 1,2$	$d = 2,3$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

Operating Instructions

CertoClav Vacuum Pro 12-22L

Laboratory autoclave

Summary

1	Validity of the operating instructions	68
2	Meaning of the labels	68
3	General Safety Instructions	69
4	Substances which must NOT be autoclaved	70
5	Introduction	71
5.1	Field of application	71
5.2	Technical data	71
5.3	Requirements for transport and storage	71
5.4	Dimensions and weight	72
6	Scope of delivery	72
7	Installation of the autoclaves	73
7.1	Packaging	73
7.2	Installation	73
7.3	Installation of accessories	76
7.3.1	Connecting the log printer	76
7.3.2	Connecting the barcode scanner	76
7.3.3	Connecting the label printer	76
7.3.4	Connecting the private CertoServer	76
7.4	Pressure sensor calibration	77
7.4.1	Calibration instructions:	77
7.5	Validation and calibration	78
8	Sterilization programs	78
9	Basic operation	79
9.1	Switching on the touch screen	79
9.2	Switching off the touch screen	80
9.3	Stand-by	80
10	Instructions for daily use	80
10.1	Preparations before use	80
10.2	Refilling of distilled water	80
10.3	Steps to start a sterilization process	81
10.3.1	Loading of the autoclave	81
10.3.2	Close autoclave door	81
10.3.3	User registration and login	82

10.3.4	Adding, Editing and Deleting programs	84
10.3.5	Select sterilization program on screen	86
10.3.6	Status screen.....	87
10.3.7	CertoTrace – Documentation and loading management	87
10.3.8	Step-by-step assistant.....	88
10.3.9	Material test mode and automatically repeated start	88
10.4	Unloading the autoclave	89
10.5	Emptying the service water tank	89
10.6	PC software.....	90
10.6.1	CertoCloud software features	90
10.6.2	Display of sterilization protocols via web browser.....	91
11	Control panel overview	92
11.1	Information menu.....	92
11.1.1	Tutorial videos	92
11.1.2	Send e-mail to CertoSupport	92
11.1.3	TeamViewer remote support and screen sharing	93
11.2	Logging menu.....	94
11.2.1	View logs	94
11.2.2	Print out protocols.....	94
11.2.3	Example print out	95
11.2.4	Tracing label printing	96
11.2.5	Tracing of marked products.....	96
12	Settings	97
12.1	User administration	97
12.2	Network settings.....	97
12.2.1	Connection to a Wi-Fi network.....	97
12.2.2	Connection to a LAN network.....	98
12.2.3	Network communication	98
12.2.4	Connection of the autoclave with the local CertoServer	99
12.3	Sterilization settings	99
12.3.1	Real-time logging on USB thumb drive.....	99
12.3.2	Step-by-step assistant.....	100
12.3.3	Export all protocols to USB thumb drive	100
12.3.4	Export all protocols to SD card	101

12.3.5	Export all protocols via network.....	101
12.3.6	Switch preheating on or off	102
12.3.7	Permanently switch preheating on or off.....	102
12.3.8	Material test mode – Automatically repeated start.....	103
12.3.9	CertoTrace documentation solution	103
12.3.10	Label designer	104
12.4	System settings	104
12.4.1	Change date and time.....	104
12.4.2	Install updates via the internet.....	105
12.4.3	Install updates via USB thumb drive.....	106
12.4.4	Install updates via SD card	106
12.4.5	Restore factory settings.....	107
12.5	Audit Trail	107
12.5.1	Audit Trail Event ID's	108
12.5.2	Change language.....	108
12.5.3	Change admin password.....	109
12.6	Configuration of SMS and e-mail notifications.....	109
12.7	GLP and GMP data	109
12.7.1	Header	109
12.7.2	Autoclave designation	109
12.7.3	Project name.....	110
12.7.4	Advanced GLP and GMP settings.....	110
12.8	Lock settings	110
13	Error messages	110
13.1	Error in sterilization process.....	110
13.2	Error in status bar	113
14	Maintenance.....	113
14.1	Adjusting door closing pressure	114
14.2	Cleaning the door seal	114
14.3	Replacing the door seal	114
14.4	Replacing the safety valve	115
15	Transport and storage	115
15.1	Preparations	115
15.2	Emptying the water tanks.....	115

15.3	Storage conditions	116
16	Warranty	116
17	Appendix	116
17.1	Appendix 1: Structure diagram	116
17.2	Appendix 2: Electronic schematic	117
17.3	Appendix 3: EMV	118

CE declaration of conformity

This is to certify that the **CERTOCLAV** autoclaves of type CV-Vac-Pro-8L, CV-Vac-Pro 12L, CV-Vac-Pro-18L, and CV-Vac-Pro-22L comply with the following EU Directives:

2014/35/EU Low Voltage Directive

Related standards:

EN 61010-2-040:2015 Safety requirements for electrical equipment for measurement, control and laboratory use
Laboratory equipment – Part 2-040: Particular requirements for sterilizers and cleaning machines
Disinfection devices for the treatment of medical material
EN 61010-1:2010 Safety requirements for electrical equipment for measurement, control and laboratory use
– Part 1: General requirements

The autoclave is designed according to the low voltage directive 2014/35/EU.

2014/30/EU Electromagnetic compatibility / EMC Directive

Related standards:

EN 61326-1:2013 Electrical equipment for measurement, control and laboratory use - EMC requirements – Part 1: General requirements
EN 61000-3-3:2013 Electromagnetic compatibility (EMC) – Part 3-3: Limits – Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply networks for devices with a rated current ≤ 16 A per conductor that are not subject to any special connection conditions

2011/65/EU Restriction of the use of certain hazardous substances (RoHS) Directive / Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment

Related standards:

EN 50581:2012 Technical documentation for the assessment of electrical and electronic equipment regarding the restriction of hazardous substances

Manufacturing quality:

CertoClav manufactures according to modern production methods using high quality materials. Decades of experience and experienced employees guarantee high quality products and maximum customer service. All devices are subject to 100% testing for functionality and quality.

Manufacturer:

CERTOCLAV Sterilizer GmbH, A-4060 Leonding, Austria

Michael Dirix
Managing Partner

Leonding, March 2020

1 Validity of the operating instructions

	This manual is only valid for CertoClav Vacuum Pro Series autoclaves with the model designations CV-08-V-Pro, CV-12-V-Pro, CV-18-V-Pro and CV-22-V-Pro.
	Please read this operating manual before the initial operation of the autoclave.
	Strictly adhere to all safety instructions in this operating manual.
	Keep these operating instructions within easy reach of the autoclave.
	If there are questions or problems with the autoclave, contact your dealer or the manufacturer.
	Operation and maintenance of the autoclave may only be carried out by trained personnel.

2 Meaning of the labels

	ATTENTION!
	PROTECTIVE GROUND WIRE CONNECTION
	Hot surface, keep away
	ENVIRONMENTAL PROTECTION – Do not dispose of electronic waste together with household waste. Contact your dealer or a local recycling center to dispose of the waste properly.
	MANUFACTURER
	COMMUNAUTÉ EUROPÉENNE
	DATE OF PRODUCTION
	SERIAL NUMBER
	THIS SIDE UP
	SHELTER FROM RAIN
	DO NOT ROLL
	STACKING: MAX. 3 UNITS
	ROOM TEMPERATURE: 5 °C~40 °C
	RELATIVE HUMIDITY ≤ 80%

3 General Safety Instructions

- ⚠ Only use the equipment in accordance with these instructions. Store this operating manual in a place that is easily accessible to all users.
- ⚠ Use a power outlet with a protective earth conductor and make sure that the protective earth conductor is well grounded.
- ⚠ Do NOT install the autoclave in a place where the power supply cannot be quickly disconnected.
- ⚠ Operate the autoclave exclusively with the supply voltage declared on the autoclave's type plate.
- ⚠ Do NOT touch the power plug or socket with wet hands.
- ⚠ Do NOT load, bend or twist the cables. Do NOT place heavy objects on the cables.
- ⚠ Do NOT autoclave pressure-tight containers; these containers could implode or explode.
- ⚠ Only fill the water tank with distilled or demineralized water.
- ⚠ Do NOT block the door of the autoclave.
- ⚠ Do NOT cover the ventilation openings at the back of the autoclave.
- ⚠ If you notice any unusual odors or noises during operation, immediately disconnect the autoclave from the power supply and contact your local dealer or the manufacturer by e-mail.
- ⚠ If you do NOT operate the autoclave for a long period of time, please disconnect the power supply and empty the internal water tanks.
- ⚠ Only connect original accessories to the autoclave.
- ⚠ Do NOT modify the device.
- ⚠ The autoclave must NOT be installed in an area, that is subject to the danger of explosions.
- ⚠ The autoclave may only be used in laboratories. It is NOT suitable for domestic use.
- ⚠ Please fill bottles and other containers only up to a maximum of half the volume with liquid.
- ⚠ Do NOT remove containers with liquids from the chamber until a liquid temperature of less than 80 °C is measured. Otherwise, there is a danger of boiling delay of the liquid.
- ⚠ If the pressure is quickly released during autoclaving liquids due to an unscheduled program termination, there is a possibility of boiling delay of the liquid. Therefore, before opening the autoclave, always wait until the temperature of the medium has reached less than 80 °C.
- ⚠ The autoclave must NOT be transported when filled.
- ⚠ The autoclave must NOT be forced open.
- ⚠ This autoclave may only be serviced by trained dealers.
- ⚠ Only use original spare parts for maintenance.

4 Substances which must NOT be autoclaved

Never use the autoclave to sterilize any of the following materials or alkaline solutions. Sterilization of such substances can lead to explosion, corrosion of the sterilization chamber and piping, as well as attrition of the seals.

Explosive substances

Nitroglycol, nitroglycerine, nitrocellulose and other explosive nitric acid esters.

Trinitrobenzene, trinitrotoluene, picric acid and other explosive nitro compounds.

Peracetic acid, methyl ethyl ketone peroxide, benzoyl peroxide and other organic peroxides.

Flammable substances

Metallic lithium, potassium, sodium, yellow phosphorus, phosphorus sulfide and red phosphorus.

Celluloids, calcium carbide (carbide), calcium phosphide and magnesium powder.

Aluminum powder, magnesium powder and metal powder other than aluminum powder.

Sodium dithionite (or sodium hydrosulfite).

Oxidizing agent

Potassium chlorate, sodium chlorate, ammonium chlorate and other chlorates.

Potassium perchlorate, sodium perchlorate, ammonium perchlorate and other perchlorates.

Potassium peroxide, sodium peroxide, barium peroxide and other inorganic peroxides.

Potassium nitrate, sodium nitrate, ammonium nitrate and other nitrates.

Sodium chlorite and other chlorites.

Calcium hypochlorite and other hypochlorites.

Combustible materials

Ethyl ether, petrol, acetaldehyde, propylene oxide, carbon disulfide and other substances whose flash points are in the range of -30 to 0 °C.

Methanol, ethanol, xylene, benzyl acetate (or amyl acetate) and other substances whose flashpoints are between 0 and 30 °C.

Kerosene, gas oil, turpentine oil, isopentyl alcohol (or isoamyl alcohol), acetic acid and other substances with a flash point between 30 °C and 65 °C

Combustible gas (hydrogen, acetylene, ethylene, methane, ethane, propane, butane and other substances operating at a temperature of 60 °C at an absolute pressure of 0.2 bar)

If a saline solution has been spilled inside the chamber, remove all liquid from the chamber and wipe the liquid around the door seal. Clean the door seal using distilled or demineralized water.

5 Introduction

This autoclave is an automatic steam sterilizer with vacuum pump and steam generator. Depending on the program and settings, the sterilization temperature is between 105 °C and 134 °C. Solid materials and liquids can be sterilized. The autoclave is outfitted with a flexible media sensor, to measure the temperature directly within the medium. The program parameters and sterilization programs are displayed on the front LCD screen and on the mobile touch screen. If a problem is detected during the sterilization process, the autoclave automatically issues an error message. The autoclave is equipped with internal water tanks for fresh and process water.

5.1 Field of application

This autoclave is exclusively intended for the use in laboratories. For example, for the sterilization of liquids, solids, waste and textiles such as swabs, glasses, petri dishes, porous materials, powders, hollow bodies and instruments. No pressure-tight sealed containers may be sterilized.

5.2 Technical data


The working area of the autoclave:

- Ambient temperature: 5 °C – 40 °C
- Relative air humidity: ≤ 80%
- Air pressure: > 70kPa;
- Supply voltage: 230V AC, 50Hz, 2500VA
- Maximum working pressure: 0.21~0.23MPa
- Maximum operating temperature: 135 °C ~ 137 °C
- Expected life: 10+ years

5.3 Requirements for transport and storage

- Ambient temperature: 5 °C ~ 40 °C
- Relative air humidity: ≤80%.
- No corrosive gases in the environment

5.4 Dimensions and weight

Model	A	B	C	D	E	F	G	H
CV-VAC-PRO-12L	415mm	503mm	433mm	540mm	610mm	200mm	128mm	178mm
CV-VAC-PRO-18L	415mm	503mm	433mm	540mm	610mm	200mm	128mm	178mm
CV-VAC-PRO-22L	415mm	503mm	433mm	640mm	710mm	200mm	128mm	178mm

Model	Chamber volume	Chamber diameter	Chamber depth	Net weight
CV-VAC-PRO-12L	12 liters	200mm	360mm	51kg
CV-VAC-PRO-18L	18 liters	250mm	350mm	53kg
CV-VAC-PRO-22L	22 liters	250mm	450mm	55kg

6 Scope of delivery


Index	Description	Quantity
1	QuickConnect discharge hose	1
2	Drawer set consisting of: 1 x handle, 3 x drawers, 1 x drawer holder	1
3	Power cable for autoclaves	1
4	Operating instructions	1
5	Mobile touch screen	1
6	Power cable for touch screen	1
7	Base for touch screen with four fastening screws	1
8	Interface box for touch screen	1
9	Autoclave	1
10	Warranty card	1
11	Data cable (1m length) between autoclave and touch screen	1

7 Installation of the autoclaves

7.1 Packaging


This autoclave is delivered in a wooden crate. To open, the metal flaps on the top of the box must be opened. Lift the autoclave vertically out of the box with the help of a second person at the two loops.

(see picture)

7.2 Installation

If you place this autoclave on a laboratory table, check that the table can support the weight of the autoclave. Please keep at least 10 cm of free space around the autoclave. Leave at least 20 cm of free space at the ventilation slots on the right side of the autoclave. Make sure that the room is well ventilated. Use a level to check that the autoclave is horizontally leveled. **Important: Place the autoclave at the front of the chamber 1.5 cm higher than the back, to achieve a slight tilt, so that the condensate in the chamber can flow towards the drain.**

Open the door of the autoclave and remove the tray set including accessories. Remove the packing material.

Mount the base to the touch screen using the four screws provided.

Now, the touch screen with the base should look like the illustration on the right.

Connect the touch screen to the interface box.

Connect the power supply unit of the touch screen to the interface box.

Only place the touch screen on flat, stable surfaces. You can also place the waterproof touch screen on top of or next to the autoclave.


Make sure that the table can support the weight of the autoclave.

Connect the power cable to the device socket at the back of the autoclave. Note: Check, whether the information on the autoclave type plate matches the information on the voltage source. Check that the safety switch of the autoclave is in the upper position. Switch the autoclave on at the green rocker switch.

Connect the autoclave with the RS232 cable at the COM 4 connection of the interface box to the RS232 connection on the back of the autoclave.

Turn on the touch screen.

If a login window appears on the screen after a few seconds, the touch screen is correctly connected to the autoclave.

If a message about a connection termination with the autoclave appears on the screen, the touch screen and autoclave are not connected correctly. Check whether the autoclave is switched on and connected to the COM4 port of the interface box.

7.3 Installation of accessories

7.3.1 Connecting the log printer

Connect the log printer (thermal or dot matrix printer) to the COM 1 connection of the interface box. The connecting cable can be found in the carton of the log printer. Follow the operating instructions of the printer.

7.3.2 Connecting the barcode scanner

Connect the barcode scanner to the USB port of the interface box. You can take the connection cable from the carton of the scanner. Follow the operating instructions for the barcode scanner.

7.3.3 Connecting the label printer

Connect the label printer to the COM2 connection of the interface box. You can take the connecting cable from the carton of the label printer. Follow the operating instructions for the label printer.

7.3.4 Connecting the private CertoServer

If you do not have/ want access to the Internet and still want to use CertoCloud, it is recommended to connect the local CertoServer to the autoclave. In chapter 12.2.4 you will learn how to connect the CertoServer to the autoclave.

7.4 Pressure sensor calibration

Since the atmospheric pressure depends on the installation location, you should recalibrate the pressure sensor after changing the installation location.

7.4.1 Calibration instructions:

Switch on the autoclave and open the door. Log in to the administrator account. User account: Admin, password: "admin". If the pressure display shows a pressure that deviates from 0kPa, the pressure sensor offset should be corrected. Navigate to Settings->Calibration and enter a suitable offset for the pressure sensor in [kPa].

Example:

The screenshot shows the 'Settings' menu with 'Calibration' selected. The 'Pressure Sensor' offset is set to -10 kPa / bar. The status bar at the bottom displays 'PRESS: 2 kPa' in red.

For example, the pressure display shows a pressure of 2 kPa when the door is open. This can be seen in red in the status bar below. Click on the value in the status bar.

The screenshot shows the 'Settings' menu with 'Calibration' selected. The 'Pressure Sensor' offset is set to -12 kPa / bar. The status bar at the bottom displays 'PRESS: 2 kPa' in red.

If the pressure sensor value in the status bar is positive. To correct the current value, enter minus the value from the current pressure. If the pressure in the status bar is in the minus range, you have to move the calibration up.
(Here 2 in plus, so the value -10 is decreased by 2: -12)

The screenshot shows the 'Settings' menu with 'Calibration' selected. The 'Pressure Sensor' offset is set to -12 kPa / bar. The status bar at the bottom displays 'PRESS: 0 kPa' in red.

Calibration was successful when 0 kPa is displayed on the bar below. Please press the button below to confirm the entries. The parameters are permanently adopted.


If the pressure sensor is not properly calibrated after the autoclave has been installed, the autoclave may refuse to start a program.

7.5 Validation and calibration

If you have any questions about validation and calibration of the instrument, please contact the manufacturer.

8 Sterilization programs

Program	Vacuum phases	Sterilization time [min]	Sterilization temperature [°C]	Sterilization pressure [kPa]	Vacuum hold time [min]	Drying time [min]
121 °C / SOLID	1	20	121	110	-	3
121 °C / POROUS	3	20	121	110	-	7
121 °C / HOLLOW	3	20	121	110	-	10
134 °C / SOLID	1	4	134	210	-	3
134 °C / POROUS	3	4	134	210	-	7
134 °C / HOLLOW	3	4	134	210	-	10
B-D TEST	3	4	134	210	-	7
VACUUM TEST	-	-	-	-80	15	-
CLEAN PROCESS	3	5	105	20	-	10
PRIONS	3	19	134	210	-	10
LIQUID	-	20	121	110	-	0


- You can create a new custom program or modify an existing program at any time.
- The SOLID programs can only sterilize solids that are not wrapped. Such as unwrapped forceps, scissors, tweezers, trays, etc.
- The "Porous" programs are specially designed for bodies made of porous material. Three deaeration steps by means of the vacuum pump ensure complete penetration of the steam into the porous material. This is a necessary condition for the proper sterilization of porous materials.
- The "Hollow" programs can sterilize type A and type B hollow items. These are, for example, glasses or cannulas.

- The "B&D/Helix" program tests hollow bodies of type A and textiles. The "B-D Test" and "Helix Test" programs are almost identical. The B&D test tests the penetration of textiles. In the Helix test, the penetration of hollow bodies of type A is tested. Thus, different indicators are used.
- The "Vacuum Test" program performs an air tightness test.
- The "Clean Program" program cleans the autoclave. If the autoclave displays the "Need Clean" message, please run this program. This program can also be used to clean other bodies.
- The "Prions" program is specially designed for the sterilization of prions on surfaces of all kinds.
- The "Liquid Program" is used for sterilizing liquids in bottles and other open containers. During sterilization, the media sensor is placed in the liquid medium or in a reference glass containing distilled water. The autoclave now adjusts the holding time during sterilization so that the medium has maintained the required sterilization temperature for at least 20 minutes. If the media sensor is not used, the media temperature cannot be logged during sterilization. CAUTION: Please never place pressure-sealed bottles, pouches or containers without a pressure equalization mechanism provided in the autoclave. These could burst due to overpressure during cooling. NOTE: Only fill jars to a maximum of half their capacity.


	Please never place tightly closed bottles, bags or containers in the autoclave. These could burst due to overpressure.
	Please perform a vacuum test every month. If the program does not complete successfully, contact your dealer or the manufacturer.
	You can store up to 50 different programs in the autoclave

9 Basic operation

Click and wipe

Long click and drag

9.1 Switching on the touch screen

1. Press and hold the -button on the front panel (to the right of the screen) for at least 3 seconds.
2. Switch on the autoclave at the green power switch.

9.2 Switching off the touch screen

Press the -button on the front panel (to the right of the screen) for at least 3 seconds. A notification will then appear on the screen. Select "Shutdown" to turn off the screen.

9.3 Stand-by

Briefly push the -button next to the screen. The backlight of the screen will then turn off. Push the -button again to reactivate the touch screen.


10 Instructions for daily use

10.1 Preparations before use

Switch on the autoclave at the green power switch. Switch on the touch screen directly on the screen. To do this, press the -button on the front of the screen (to the right of the screen) for at least 3 seconds. The autoclave's touch screen then switches on and displays the login window. The autoclave is ready for operation.

	Only use the device with the mains voltage specified on the type plate.
	Before the first use, the touch screen displays the message "Please fill distilled water into the water tank". As soon as the tank is filled, the notification disappears, and the autoclave is ready for operation.

10.2 Refilling of distilled water

When the message "Please fill distilled water into the water tank" appears on the screen, fill distilled water into the water tank at the top of the autoclave. As soon as the warning message disappears and a beep sounds, the minimum fill level of distilled water has been reached. Take care not to overfill the tank. The inlet to the water tank is located on the lid of the autoclave. Only distilled or demineralized water is suitable for the operation of the autoclave.

If the tank is overfilled, water can be drained with the QuickConnect hose connection at the front of the unit. If the message "Please fill distilled water into the water tank" appears on the screen during a sterilization run, you do not need to cancel the sterilization run. However, please fill distilled water into the tank before starting the next sterilization cycle.

	Neither the manufacturer nor your dealer is liable for damage caused by using contaminated water.
	Only fill demineralized or distilled water into the freshwater tank. A water quality sensor continuously monitors the quality of the water. If too impure water is filled into the tank, a warning message appears on the front integrated LCD screen and the touch screen.

10.3 Steps to start a sterilization process.

Before starting the program, first check whether the autoclave is ready for operation. If no warning or error message is displayed on the touch screen, the autoclave is ready for operation.

10.3.1 Loading of the autoclave

Please make sure that the total volume of the load is less than 70% of the volume of the sterilization chamber. If you put your load in or take it out of the sterilization chamber, please use the tool provided.

	Please leave at least 10 mm free space between the load and the sterilization chamber wall.
	When running the B&D/Helix program, place the indicator in the inner center of the sterilization chamber. This way you get a reliable test result.

10.3.2 Close autoclave door

	When the temperature in the sterilization chamber is high, the autoclave door is difficult to close. Wait 30 seconds before closing the door.
---	---

10.3.3 User registration and login

After switching on the touch screen, the login window appears. Please select an user within the spinner. After you insert the right password you can proceed with the login by pressing the Login-Button below.

STEAM: 79 °C MEDIA: 76 °C PRESS: 1 kPa support@certoclav.com

To create a new user, first log in with a user with admin rights and navigate to the settings area. Press the plus symbol at the top right of the screen. This becomes visible as soon as you are in the "User Account" tab.

A dialog appears. To create a local account without a CertoCloud server connection, press the left button of the dialog. To add an existing CertoCloud account to the device, press the Existing button. If you don't have a CertoCloud account yet, press the Create New button. To create a new CertoCloud account, an internet connection via ethernet or Wifi is required.

To create a local account (without server connection), enter your user data. The text fields can be selected and edited by clicking on the text field. Please choose a password with at least four characters and retype it. Confirm the user registration by clicking the Register User button. To abort the registration, press the Back button. To complete the registration, press the Register button

To add an existing CertoCloud account, enter your e-mail address and password. Afterwards press Add Account to confirm the account. To abort the registration, press the Back button on the left-top of the screen.

To create a new CertoCloud account, enter your user data into the appropriate text fields. Please enter your phone number with the international prefix. To validate your phone number and e-mail address, an e-mail and an SMS will be sent to your e-mail and your phone. Please enter your password twice and press the Register button.

To login with the created user account, choose your email address from the list and enter your password. Please make sure that the online mode is enabled. The online mode settings can be found under: Settings->Network->Online Mode.


If you log in for the first time, you will be asked for an activation key. You received the activation key via e-mail and via SMS, after registering for the CertoCloud account. Please enter the activation key in the text field and press the Activate Account button. Your CertoCloud account is now activated and you can use the CertoCloud account with your autoclave.

10.3.4 Adding, Editing and Deleting programs


The Certo Clav Vacuum Pro Series contains 11 ready-made programs (as described in 10.3.1). These can be changed and deleted. In addition, there are specially defined programs whose creation is dealt with below. The desired name, the sterilization time and temperature, the drying time and the vacuum steps can be freely determined.

Add Program

When you click the button on the top right with the “+” symbol, a window opens for entering the required data.


You can change this using the plus and minus buttons to the level of the value to be changed. You can choose between different vacuum units in the range from 1 to 10. The sterilization temperature can also be modified specifically. Temperatures from 105 to 134 degrees Celsius are available for selection. The sterilization times can be selected for you in the range of 4-60 minutes. The drying interval can last a maximum of 25 minutes.


By clicking the text field for the program name at the top left, you can specify any name using the input keyboard.

For example, if you want to create a program for a liquid, you have to select the selection menu on the top right. The setting options that are not required are then removed.


In order to then save a program in the database and add it to the selection menu, press the "Save Program" button. You can cancel the process at any time using the back button at the top left.

Edit Program


Of course, in addition to creating new programs, you can also edit existing programs. This refers both to the programs which are delivered in the factory state as well as to self-defined programs of the autoclave. When you click on the pencil symbol in the top right of a program window, a new window with editable program parameters opens.


Delete Program

Both factory programs and user-defined programs can be deleted at any time. To do this, go to the pencil symbol of a program of your choice.



A window then opens with a "Delete Program" button on the bottom left of the page.


10.3.5 Select sterilization program on screen

1. After logging in (see section 10.4.) navigate to the sterilization menu.
2. In the middle 4 programs are visible. Scroll down the list to see all available programs.
3. Press one program to select it. Confirm that you want to start the program.
4. If enabled in the settings, step-by-step videos are displayed before the program starts.

10.3.6 Status screen


During the current cycle, the status screen is visible. On the left side, information about the current cycle is listed. On the right side of the screen, the current steam temperature and pressure are displayed. If the liquid program has been started, the liquid temperature is also visible. To display a graphical view of the temperature and pressure curve, hold and drag the autoclave icon to the left. Afterwards a graphical recording becomes visible.

10.3.7 CertoTrace – Documentation and loading management

CertoTrace is a software that enables simple and professional documentation of the entire sterilization process and traceability of sterilized goods.

- 1) If CertoTrace is enabled (Settings->Sterilization->CertoTrace Documentation), the CertoTrace wizard is executed immediately after a program is selected.
- 2) Scan the barcode of bottles, instruments or other items placed in the autoclave to trace their sterilization with CertoTrace.


- 3) After an item has been scanned, it appears in a list in the CertoTrace scan menu. The object and the barcode are recorded in the protocol of the sterilization process.
- 4) After all items have been scanned and placed in the sterilization chamber, press the "Next"- button.
- 5) To cancel the CertoTrace wizard, press the "Back"- button in the upper left corner of the screen.
- 6) If no barcode exists for an item, use the label printer designer to create and print a barcode for the item. The label printer designer can be found under Settings->Sterilization->Label Printer Designer.

10.3.8 Step-by-step assistant


To activate the step-by-step assistant, navigate to the settings menu by clicking on the  icon in the upper right corner of the screen. Select the Sterilization option on the left side of the settings window. Press Step-by-Step Assistant to enable or disable it.

The step-by-step assistant guides the user through all necessary steps to operate the autoclave correctly by means of video instructions.

10.3.9 Material test mode and automatically repeated start

The material test mode can be activated via the checkbox in Settings->Sterilization->Material test. If the material test mode has been activated, the autoclave asks, before the program starts, how often the program should automatically repeat cycles.

A number between 1 and 1000 can be entered. The program will execute the set number of cycles fully automatically, one after the other. The status of the currently running program can be called up at any time, either browser-based at www.certocloud.com or offline via the CertoServer at <http://certocloud>. If you have activated e-mail and SMS notifications, you can retrieve data on the current status of the material test via SMS or e-mail. The e-mail and SMS notifications can be configured in Settings->Notifications. However, SMS notifications are only possible with the CertoCloud Premium account.

It is necessary to fill the freshwater tank regularly and empty the waste water tank regularly.

10.4 Unloading the autoclave

Once the sterilization cycle is complete, the screen will display "Program Complete" and you will hear three beeps. Now you can open the door and take out the sterilized load.

	Please use the supplied accessories to remove the load from the sterilization chamber. Do not open the door until the pressure inside the sterilization chamber is in the range of -5kPa to 5kPa.
	To independently check the effectiveness of the sterilization cycles, place indicators together with the load in the sterilization chamber.

10.5 Emptying the service water tank

If the screen displays "Empty the hot water tank", the wastewater tank is full. Please empty the process water with the QuickConnect hose supplied.

	Caution, the process water may be hot.
---	--

10.6 PC software

10.6.1 CertoCloud software features

The CertoCloud is a browser-based management software for your autoclave. It is available for PC, smartphone and tablet at www.certocloud.com. The CertoCloud is also available as an in-lab private cloud. For this purpose, a CertoServer must be connected to your laboratory network. The CertoServer does not require an Internet connection.

Create a user account, directly at the autoclave or web-based.


Tracking of real-time data from the autoclave.

Management of SMS and e-mail notifications. For example, when maintenance is due. *

Download all sterilization protocols as signed PDF, Excel and ZIP files.

Cycle	Program Name	Start
76	cloud generated	7/22/ PM
75	121°C SOLID	7/22/ PM
74	Neues Programm	7/22/ PM
74	121°C SOLID	7/22/ PM

View all logs in the graphical view, as a list or in detailed view.


Prove the validity of your exported protocols.

*) SMS notifications are only possible with the CertoCloud Premium account!

10.6.2 Display of sterilization protocols via web browser


If your autoclave is connected to the Internet or to the local CertoServer, you can view and manage your protocols. A sample protocol is shown below.

11 Control panel overview


11.1 Information menu

Navigate to the “Information”-tab to display the Information menu. Here, it is possible to watch training videos, contact CertoSupport directly from the touch screen, or order spare parts, accessories or services for the autoclave.

11.1.1 Tutorial videos


To scroll through the list of all available videos, press the “Video” section of the “Information”-tab. Select the desired video, to view it in full screen mode.

(Example picture)

11.1.2 Send e-mail to CertoSupport


To send an e-mail to CertoSupport directly from the device, the device must be connected to the Internet. To connect the device to the Internet, you can plug a network cable into the Ethernet port or establish a wireless connection to your wireless router. Press the “SUPPORT”-button to open the e-mail window.



Select the subject of your message to CertoSupport from the drop-down menu. Then enter your message to CertoSupport in the message text field. To send the e-mail, press the “Send”-button.

11.1.3 TeamViewer remote support and screen sharing


To split the screen with CertoSupport using TeamViewer, navigate to the “Information”-tab and press the “Support”-button. Remote support is only possible by making an appointment by e-mail or telephone beforehand.

As soon as you start a TeamViewer session, a window appears which gives you the opportunity to roughly describe the problem once more. Then press “Start”.

After you have pressed “Start”, wait until CertoSupport has joined the TeamViewer Session. Once CertoSupport has joined the session, they will ask for your permission to share the screen. Press “Allow”. Now the technician can access the menu of your device. To cancel the screen share, drag the “TeamViewer” icon to the “X” icon at the bottom of the screen.


TimeStamp [min]	Chamber temp. [°C]	Media temp. [°C]	Pressure [kPa]
1.00	48.92	0.0	2.00
2.00	48.59	0.0	2.00
3.00	47.89	0.0	2.00
4.00	47.89	0.0	2.00
5.00	47.52	0.0	2.00
6.00	47.14	0.0	2.00
7.00	46.76	0.0	2.00
8.00	46.36	0.0	2.00

11.2 Logging menu


11.2.1 View logs

To view sterilization protocols, navigate to the “Logging”-tab. A list of all recorded protocols is displayed on the left side of the screen. The “Cloud” icon in the protocol list shows whether the protocols have been backed up to the CertoCloud or to your local CertoServer. Press the drop-down menu above the protocol list to sort the protocols by specific criteria. To view the details and a graphical view of a protocol, select it by touching the respective list entry. On the right side of the window, the graphical history of the protocol is now displayed. To display the protocol as a list, press the “List” icon above the graph.

11.2.2 Print out protocols

- 1) Make sure that a log printer is connected to the device.
- 2) Navigate to the “Logging”-tab.
- 3) On the left side of the screen, select the protocol you want to print.
- 4) Press the „Printer” icon in the upper right corner of the screen.
- 5) Select the “Print protocol” option.

11.2.3 Example print out

Optional


```

Header line 1
Header line 2
Header line 3

Protocol      CertoClav Vacuum Pro Series
Model:        VAC-PRO-22L
S/N:          22ABCDEFGHJ
Program:       134°C HOLLOW
Program description:
    3 x venting
    134 °C    2.1bar    30min
    Hollow instruments
    wrapped or unwrapped
Media:         NO
Cycle no:      971
Start time:    12.01.15 17:42:00
End time:      12.01.15 18:16:00
Result:        PASSED
hh:mm:ss      Temp °C    Press [bar]
17:42:00      27.5       1.00
17:42:20      27.5       0.50
17:42:40      27.6       0.18
17:43:00      80.4       0.80
[...]
18:15:00      65.9       0.18
18:15:20      65.4       0.18
18:15:40      68.0       0.50
18:16:00      70.2       1.00
Program finished successfully
END

Signed by _____
Confirmed by _____

```


11.2.4 Tracing label printing

- 1) Make sure a label printer is connected to the machine.
- 2) Navigate to the "Logging"-tab.
- 3) Select the desired protocol on the left side of the screen.
- 4) Press the „Printer" icon in the upper right corner of the screen.
- 5) Select the "Print a tracing label" option.
- 6) Enter the shelf life of the sterilized products.
- 7) Press "OK" to print the self-adhesive label with expiry date and tracing barcode.

11.2.5 Tracing of marked products

- 1) Navigate to the Logging tab.
- 2) Press the "Scanner" icon, next to the "Printer" icon.
- 3) A dialog box appears. Now scan the barcode of the sterilized product using the barcode scanner.
- 4) The protocol of the sterilized product now appears on the screen. On the right side of the screen a graphical and list view of the protocol can be viewed.

12 Settings


Press the  icon at the top right of the screen to enter the “Settings” menu.

The following settings can be adjusted here:

12.1 User administration

All registered users are now displayed on the right side of the window. Here you can change the user data and passwords via the “Edit”-button and remove a user account from the touch PC by pressing the “Delete”-button.

Add users by clicking the button in the upper right corner.


See chapter 10.3.3 for detailed instructions on user registration and login process.

12.2 Network settings


12.2.1 Connection to a Wi-Fi network

1. Click on the “Network”-button.
2. If the Wi-Fi checkbox is blank, select the empty checkbox to enable Wi-Fi.
3. Then press “Manage WLAN connections” to display a list of available WLAN networks. Select the relevant WLAN network and enter access keys and other parameters, if necessary.

12.2.2 Connection to a LAN network


The autoclave can also be connected to a LAN network via ethernet cable. Connect the network cable to the ethernet port of the interface box.

12.2.3 Network communication



If the autoclave data is to be synchronized with the CertoCloud (www.certocloud.com) or the local CertoServer (<http://certocloud>), the online mode must be activated.

1. Select the "Network"-button on the left side of the settings window.
2. Click the "Online Mode" option to enable or disable data synchronization with CertoCloud. **Note:** If the online mode is deactivated, your protocols are only stored locally on the autoclave. As soon as you reactivate online mode, all protocols stored locally up to that point, are uploaded to CertoCloud and are available online in your backup storage.



12.2.4 Connection of the autoclave with the local CertoServer

1. Connect the CertoServer to the interface box or the in-house network.
2. Connect the interface box to the same network as the CertoServer.
3. Navigate to: Settings->Network->Select server.
4. Select the option "Add server from local network".

5. The autoclave searches for the CertoServer in the local network.

6. After the CertoServer has been found, confirm with the "Save"-button.

12.3 Sterilization settings


12.3.1 Real-time logging on USB thumb drive

You can use a USB thumb drive to save the sterilization protocols in real time during sterilization. Please connect a USB thumb drive to the front USB port of the autoclave.

12.3.2 Step-by-step assistant





1. Enable or disable the step-by-step assistant by clicking the checkbox

12.3.3 Export all protocols to USB thumb drive

The touch PC can store around 1000 protocols in its internal memory. To export all protocols to a USB thumb drive, insert a USB thumb drive into the interface box.




Navigate to Settings->Sterilization->Export protocols to USB thumb drive.

If the protocol export fails, your USB thumb drive is not being recognized immediately. In this case, leave the USB thumb drive plugged into the device and restart the touch screen. Now try to export the protocols again.


12.3.4 Export all protocols to SD card

1. Insert a micro-SD card into the hinged slot on the right side of the touch screen.
2. Click "Export all protocols to micro-SD card".
3. If the export of the protocols fails, restart the touch screen with the SD card inserted. Now try to export the protocols again.


12.3.5 Export all protocols via network

2. Log into the autoclave with your CertoCloud or CertoServer account. Make sure that the online mode is activated. If the login is successful, all logs are uploaded to CertoCloud. The export was successful, if the cloud icon is displayed with a check mark.


12.3.6 Switch preheating on or off

It is recommended to heat the chamber to 50 °C prior to the vacuum phase of the sterilization program. To activate this function, navigate to Settings->Sterilization->Chamber “Preheat ON/OFF” and activate the check box.

12.3.7 Permanently switch preheating on or off



You can also keep the sterilization chamber in stand-by mode at 50 °C operating temperature for a maximum of 6 hours, to optimize the batch time. To enable this function, navigate to Settings->Sterilization->Hold temperature ON/OFF and select the checkbox. After the maximum time has expired, the heating switches off automatically. To reactivate the function, you must restart the device. The setting remains stored in the menu.

12.3.8 Material test mode – Automatically repeated start

To enable or disable the material test mode, navigate to the checkbox under Settings->Sterilization->Material Test Mode.



If the material test mode is activated, before starting a program for the first time, you will be asked how often it should be executed. This has the advantage that longer material tests can be performed fully automatically.

12.3.9 CertoTrace documentation solution


CertoTrace is the software for simple and professional documentation of the entire sterilization process. CertoTrace records which products were sterilized, by whom and when. Traceability of sterilized products and management of expiry dates is already integrated into the workflow. The CertoClav label printer prints expiry dates and a barcode for product traceability on self-adhesive labels.

To activate CertoTrace, navigate to Settings->Sterilization->CertoTrace load documentation. Check the box to activate CertoTrace. After activation, CertoTrace is automatically integrated into the workflow when a program is started.

12.3.10 Label designer


To mark items with barcodes, you can use the integrated label designer. The label designer can be opened under Settings->Sterilization->Label Designer.


12.4 System settings

12.4.1 Change date and time


To change the date or time, navigate to the Settings menu. This is marked with  symbol at the top right of the window.


The sub-item "Date and Time" can be found under the section "Device". The current time is displayed here. When you click on this field, a dialog box for selecting the date and time appears.

On the left side you can choose whether you want to change the day or year. You can then confirm this selection by selecting the desired date with the "OK" button.


Another dialog window opens for entering the desired time. After the hour has been selected, the minute display appears automatically, which can also be edited.

You can cancel the entry by touching any point outside the dialog window or by pressing the "Cancel" button. You can confirm the entry with the "OK" button on the lower side of the dialog window.


12.4.2 Install updates via the internet

1. Make sure the autoclave is connected to the Internet.
2. Click on the "Device" button.
3. In the details view on the right side, select "Check for updates".
4. A progress bar appears with the available download. Once all updates have been downloaded, they are installed automatically. If a new software version is available, an installation wizard will guide you through the installation process.


12.4.3 Install updates via USB thumb drive

1. Download the latest update from www.certoclav.com or ask for the latest update for your device by sending an e-mail to support@certoclav.com. Please inform the support team of the serial number of the autoclave. The serial number can be found on the type plate of the autoclave or under Settings->Device->Serial number. Copy the "update.zip" to the root directory of a USB thumb drive and insert this USB thumb drive into the interface box.
2. Click on the "Device" button.
3. In the detailed view on the right-hand side, scroll down and select the option "Installation of updates via USB thumb drive".
4. If the USB thumb drive is not recognized, restart the touch screen and try again.
5. Follow the instructions on the screen to install the respective updates.

12.4.4 Install updates via SD card

1. Download the latest update from www.certoclav.com or ask for the latest update for your device by sending an e-mail to support@certoclav.com. Please inform the support team of the serial number of the autoclave. The serial number can be found on the type plate of the autoclave or under Settings->Device->Serial number. Copy the "update.zip" to the root directory of a micro-SD card and insert this SD card into the slot on the right side of the touch screen.
2. Click on the "Device" button.
3. In the detailed view on the right-hand side, scroll down and select "Install updates via SD card".
4. If the SD card is not recognized, restart the touchscreen and try again.
5. Follow the instructions on the screen to install the respective updates.


12.5.1 Audit Trail Event ID's

Event ID	Description
1	ACTION_OPEN_WINDOW
2	ACTION_USER_CREATED
3	ACTION_USER_CHANGED
4	ACTION_USER_DELETED
5	ACTION_PROGRAM_STARTED
6	ACTION_PROGRAM_INDICATOR_CHANGED
7	ACTION_PROGRAM_CANCELED
8	ACTION_PROGRAM_FINISHED
9	ACTION_PROGRAM_CREATED
10	ACTION_PROGRAM_CREATED_FAILED
11	ACTION_PROGRAM_DELETED
12	ACTION_PROGRAM_DELETION_FAILED
13	ACTION_PROGRAM_SAVED
14	ACTION_PROGRAM_SAVE_FAILED
15	ACTION_PROGRAM_FAILED
16	ACTION_SUCCESS_LOGIN
17	ACTION_FAILED_LOGIN
18	ACTION_LOGOUT
19	ACTION_ACTIVATED
20	ACTION_PREFERENCE_SETTINGS_CHANGED
21	ACTION_ERROR
22	ACTION_PREFERENCE_LOCKOUT_MENU_CHANGED

12.5.2 Change language

1. In the settings menu, open the language menu by clicking on the “Language” button.
2. Select “Language” in the detailed view on the right side of the screen.
3. You will see a scrollable list of supported languages. Select a language by clicking on the desired language.

12.5.3 Change admin password

If you forget the admin password, you will not be able to access some functions. Please choose it carefully.

1. Click on the "Device" entry
2. scroll in the detail view on the right side to the entry "Admin password" and select it.
3. A window opens, here you can set and confirm the new password.

12.6 Configuration of SMS and e-mail notifications

If you do not have a CertoCloud account, you can create a new CertoCloud account via the login window. Note: To create a CertoCloud account, the autoclave must be connected to the Internet. Log in with your CertoCloud account.

Navigate to  Settings->Notifications to configure SMS and e-mail notifications.

If you would like to receive an e-mail when a program has been terminated, maintenance is required, or a program has been successfully completed, then activate the corresponding e-mail checkbox. Enter the e-mail address of your choice, which should receive the notification. If you have a CertoCloud Premium account, you can also enter a mobile phone number to receive an SMS notification. Notifications via SMS are only possible with the CertoCloud Premium account. Contact your dealer, if you are interested in an upgrade.

Finally, click "Save" to save the changes.

12.7 GLP and GMP data

Navigate to  Settings->GLP and GMP Data to view the GLP and GMP settings.

12.7.1 Header

If this setting is activated, the header is displayed and printed at the head of each sterilization protocol. Up to four lines of text are available.

12.7.2 Autoclave designation

Set a name for the device in this setting. An alphanumeric text can be entered. The default name is empty.

12.7.3 Project name

Click the "Project Name" button to enter a project name. An alphanumeric text of up to 25 characters can be entered. The default name is empty.

12.7.4 Advanced GLP and GMP settings

Select the following checkboxes in the advanced GLP and GMP settings to configure the printout and digital copy of the sterilization protocols.

Setting	Option (default value is marked bold)
Print Autoclave Designation	YES , NO
Print User Identification	YES , NO
Print Project Name	YES , NO
Print Application Name	YES , NO
Print Signature Fields	YES , NO
Print Header	YES , NO
Print Date	YES , NO
Print Autoclave Serial Number	YES , NO

12.8 Lock settings

To prevent unauthorized changes to the autoclaves, sensor parameters, programs and settings, the settings can only be changed in the password-protected administrator user account. Here you can define whether other user roles have access to certain settings.

13 Error messages

13.1 Error in sterilization process

If an error occurs during a program cycle, the autoclave sounds an acoustic alarm and displays an error code on the screen. In addition, the autoclave automatically releases the pressure from the chamber and automatically switches off the heating. This security system ensures the safety of the users. If an alarm occurs, please write down the error code. Do not open the door until the autoclave has cooled down and is depressurized. Only then can the chamber be opened.

	We recommend starting the program again, to see if the error occurs repeatedly.
---	---

Each of the following error messages is accompanied by a long, recurrent beep.

Code	Status	Possible cause	Troubleshooting
E -1	The program was terminated manually by the user.	User pressed the "STOP" button during sterilization.	-
E -2	The communication between external touch screen and mainboard is interrupted.	The serial cable (COM 4) between touch screen and autoclave has loosened.	Check the cable connection between autoclave and touch PC. Otherwise restart the device and the touch PC. If the problem persists, navigate to Settings->Device-> Factory Reset to reset the touch PC to the factory settings.

Code	Status	Possible cause	Troubleshooting
E +2	The sterilization pressure was exceeded by more than 40kPa.	The vacuum pump is defective.	Check the vacuum pump and then perform a vacuum test.
E31	The temperature of the chamber is higher than 150 °C.	The steam generator sensor is defective.	Check, whether the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E32	The temperature of the chamber heating is higher than 280 °C.	The steam generator sensor is defective.	Check, whether the temperature sensor (PT1000) of the chamber heater is properly inserted. If necessary, replace the temperature sensor.
E51	The temperature in the chamber is less than 0 °C.	The steam generator sensor is defective or the temperature in the chamber is too low.	Check, whether the chamber temperature is above 0 °C. Check that the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E52	The temperature of the chamber heating is less than 0 °C.	The steam generator sensor is defective or the temperature in the chamber is too low.	Check, whether the chamber temperature is above 0 °C. Check that the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E63	The temperature of the steam generator is less than 0 °C. The temperature of the steam generator is higher than 230 °C.	The steam generator temperature sensor is defective. Steam generator temperature is higher than 230 °C. Steam temperature control is unstable, causing steam generator temperature to become hotter than 230 °C.	Check, whether the steam generator temperature sensor is properly installed. Check steam generator temperature sensor (PT1000). If necessary, replace the temperature sensor.
E61	The temperature control is unstable. The control temperature inside the chamber was exceeded by 6 °C.	The chamber temperature sensor is defective.	Check, whether the temperature sensor (PT1000) is properly installed. If necessary, replace the temperature sensor.
E62	The temperature of the chamber heating is higher than 155 °C. The temperature control is unstable.	The chamber temperature sensor is damaged. The mainboard is damaged.	Check, whether the chamber temperature sensor is properly inserted. Check the temperature sensor (PT1000). Check the mainboard.

Code	Status	Possible cause	Troubleshooting
E41	During the preheating period, a chamber temperature of 100 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E42	During the preheating period, a chamber temperature of 110 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E5	During the venting phase, the pressure in the chamber is still greater than 0.5 bar after 10 minutes.	Solenoid valve is blocked.	Check the water drain valve. Clean or replace if necessary.
E6	The door was opened during operation or the door sensor is damaged.	Door sensor is damaged.	Check the door sensor and micro-switch.
E7	The measured absolute atmospheric pressure is less than 70kPa.	Incorrect air pressure. The measured absolute atmospheric air pressure is too low for this location and must be corrected.	Pressure sensor calibration. See chapter 7.5.
E8	In the pre-vacuum phase, the temperature in the chamber rises by less than 3 °C within 5 minutes.	No water in the tank or air escapes from the chamber.	Check the water level in the fish water tank. Check the chamber heater and the mainboard.
E9	During the holding phase during sterilization, the measured pressure deviates by 0.3 bar from the intended pressure.	No water in the tank. Steam generator heating defective or water pump defective.	Check the steam generator heating, water pump and water level.
E10	The electronic interlock has the wrong end position.	Electronic interlock defective. Mainboard defective.	Check the electronic interlock and the mainboard.
E11	The electronic interlock has the wrong end position.	Electronic interlock defective. Mainboard defective.	Check the electronic interlock and the mainboard.

Code	Status	Possible cause	Troubleshooting
E12	The vacuum in the sterilization chamber does not reach less than -70kPa twice in a row with a program that has at least 3 vacuum phases.	Vacuum pump defective, fan defective.	Check the vacuum pump and fan.
E99	The communication between front display and mainboard is interrupted.	Ribbon cable not connected. CPU IC has loosened.	Check the flat ribbon cable and the connection plugs between front board and main board. Check whether the microcontroller IC on the mainboard is still firmly in the socket.

13.2 Error in status bar

If the pressure display in the status bar is shown in red letters, you must recalibrate the pressure sensor of your autoclave. To do this, please follow the instructions in Chapter 7.4.

14 Maintenance

Make sure the autoclave is disconnected from the power supply during maintenance and repair. Make sure the autoclave is at room temperature.

The following parts must be checked or replaced regularly:

- Bacteria filter
- Sealing ring
- Safety valve

Component	Verification by	Cycle	Maintenance
Door	Trained staff	2 years	See service instructions or chapter <i>Adjusting door contact pressure</i> (14.1.)
Sealing ring	User	1 year	See service instructions or chapters <i>Cleaning/Replacing the door seal</i> (14.2. & 14.3.)
Printer paper	User	As required	See the operating instructions of the printer.
Safety valve	Trained staff	1 year	See service instructions or chapter <i>Replacing the safety valve</i> (14.4.)

14.1 Adjusting door closing pressure

1. Pull the rotation lock lever down and hold the lever in this position (see figure).
2. Turn the round door counterclockwise to reduce the distance between the sterilization chamber and the door. As a result, the door is closer to the chamber and the contact pressure of the door is increased.
3. Turn the round door clockwise to reduce the tightness between the door and the sterilization chamber. Consequently, the door can be closed with less force.
4. Activate the rotation lock by continuing to rotate the door until it locks by itself.


	Start a vacuum test after each door adjustment. If the test is negative, the door should be readjusted.
	If the door is too loose, it is possible that the door seal will be pushed out of its socket with a loud bang. Please make sure that the door is not too loose by performing a test run with the SOLID 134 °C program.

14.2 Cleaning the door seal

The door seal should be cleaned regularly. Please use distilled water to clean the sealing ring. If water leaks still occur after cleaning, the door contact pressure must be adjusted, or the door seal must be replaced.

14.3 Replacing the door seal

1. Hold the ring seal with one hand and use your other hand to carefully lever the ring seal out of its duct with a screwdriver. Then slowly pull on the seal to remove it completely.
2. After removing the ring seal, clean the door with distilled water.
3. Before installation, also clean the new seal with distilled water and reattach it to the door.

14.4 Replacing the safety valve

The safety valve must only be replaced by trained personnel.

1. Check the safety valve regularly for visual defects. If a flaw or defect is found, the safety valve must be replaced.
2. Remove the pipe clamp (part 1 in the illustration). Then remove the hose to which the safety valve is connected.
3. Remove the nut on the hose attachment (part 2 in the illustration).
4. Remove the safety valve and replace it with a new one.

The screenshot shows the control panel with tabs for INFORMATION, STERILISATION, and LOGGING. The STERILISATION tab is active, displaying a list of cycles on the left and a detailed log on the right.

Timestamp [min]	Chamber temp. [°C]	Media temp. [°C]	Pressure [kPa]
1.00	48.92	0.0	2.00
2.00	48.59	0.0	2.00
3.00	47.89	0.0	2.00
4.00	47.52	0.0	2.00
5.00	47.14	0.0	2.00
6.00	46.76	0.0	2.00
7.00	46.36	0.0	2.00
8.00			

On the left, the cycle list shows:

- 121 °C SOLID** (ID: 27 | 21:06 | 2017-05-12 | support@certoclav.com)
- 121 °C SOLID** (ID: 26 | 20:29 | 2017-05-12 | support@certoclav.com)
- USER DEFINED** (ID: 25 | 20:24 | 2017-05-12 | support@certoclav.com)
- 121 °C HOLLOW** (ID: 24 | 17:59 | 2017-05-12 | support@certoclav.com)

The log for the current cycle (ID: 26) shows a status of "Successfully completed" at 3.00 minutes.

Please use only original spare parts from the manufacturer. Please contact your dealer if you have questions regarding maintenance and repair.

15 Transport and storage

15.1 Preparations

Allow the autoclave to cool down and remove the power plug from the socket.

15.2 Emptying the water tanks

Please empty all internal water tanks by connecting the QuickConnect hose to the fresh water drain and to the wastewater drain and drain the water.

15.3 Storage conditions

The autoclave should be stored at a room temperature between 5 °C and 40 °C maximum. The relative humidity must not be higher than 80%. Make sure that the storage area is well ventilated and does not contain any corrosive gases.

Do not drop the autoclave during transport.


16 Warranty

Guarantee: **2 years**. The entitlement to free service is forfeited if any of the following situations arise:


- The damage was caused by incorrect installation and/or operation.
- The damage was caused by negligent behavior.
- The damage was caused by force majeure such as inadmissible supply voltage or fire.

17 Appendix

17.1 Appendix 1: Structure diagram

17.2 Appendix 2: Electronic schematic

17.3 Appendix 3: EMV

Electromagnetische Emissionen

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Steam sterilizer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Vacuum	The Steam sterilizer is suitable for use in all establishments, including domestic establishments and those directly connected to public low-voltage power supply networks that supply buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	2 kV for power supply lines 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	1 kV line(s) to line(s) 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations in power supply input lines	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Steam sterilizer requires continued operation during power

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
IEC 61000-4-11	(60% dip in <i>UT</i>) for 5 cycles 70% <i>UT</i> (30% dip in <i>UT</i>) for 25 cycles <5% <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	(60% dip in <i>UT</i>) for 5 cycles 70% <i>UT</i> (30% dip in <i>UT</i>) for 25 cycles <5% <i>UT</i> (>95% dip in <i>UT</i>) for 5 sec	mains interruptions, it is recommended that the Steam sterilizer be powered by an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 150 kHz to 80 MHz 3 Vrms 3 V/m 80 MHz to 2.5 GHz 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Steam sterilizer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2,5 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level for each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Steam sterilizer is used exceeds the applicable RF compliance level above, the Steam sterilizer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Steam sterilizer. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between

Portable and mobile RF communications equipment and the Steam sterilizer

The Steam sterilizer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Steam sterilizer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Steam sterilizer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
		m	
W	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2$	$d = 1.2$	$d = 2.3$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23