
Laboratory Equipment

in customized premium quality

MADE IN GERMANY

LET'S SHAKE IT

Edmund Bühler GmbH

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Edmund Bühler GmbH

Founded by Edmund Bühler in the university city of Tübingen in 1880, Edmund Bühler GmbH is still owned by the family today.

More than 140 years later, Edmund Bühler GmbH has developed to become a company which operates throughout the world. The production and sales of conventional laboratory devices is traditionally a main focus of the company. A second focal point is the development and production of research equipment, in particular for the quickly growing field of materials science.

In the laboratories of the world, the „Bühler Shaker“ is considered to be a synonym for the efficiency of the company. Competence, experience and the commitment to excellence in quality, which is steadily anchored in the company philosophy, form the basis for a unique program of technically mature products.

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


Description

Stainless Steel Shaker especially for incubators with CO₂ concentration

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Basic equipment

Basic device incl. basic platform with rubber mat and external control unit

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Technical Data

Motion: orbital
Shaking platform: 200 x 295 mm
Max. load: 2 kg
Shaking speed: 5 - 220 rpm, in steps of 5
Stroke: 20 mm
Runtime: 0 - 99 h 59 min / continuous
CO₂ concentration: 20 %

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Electrical supply: 100 - 240 V, 50 / 60 Hz
Enclosure protection: IP 21
Power consumption: max. 4 W
Heat emission: approx. 0.5 - 1.5 W
Ambient temperature: 4°C to 45°C
Relative humidity: ~ 95 % (non-condensing)
Length of standard cable (KM CO₂): 200 cm
Length of ribbon cable (KM CO₂-FL): 120 cm
Dimensions (w x d x h): 230 x 340 x 145 mm
Weight: 7.6 kg

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Loading Capacity (Examples)

Erlenmeyer flasks 100 ml 9 pieces
Erlenmeyer flasks 500 ml 2 pieces
Further loading possibilities see page 47/48

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The external control unit with LED display enables easy operation outside the incubator. Cable feed-through required in the incubator.

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The external control unit with LED display and ribbon cable enables easy operation outside the incubator (KM CO₂-FL).

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Universal tray KM Mini
No. 0052 119

Coated tray with drillings for individual loading with spring clamps or test tube racks.
Shaking platform 200 x 295 mm.

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Universal tray KM
No. 0052 030

Coated tray with drillings for individual loading with spring clamps or test tube racks.
Shaking platform 356 x 300 mm.

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Description

Highly reliable shaker which requires little space

Basic equipment

Basic device incl. basic platform and rubber mat

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Technical Data

Motion: orbital
 Shaking platform: 200 x 295 mm
 Max. load: 2 kg
 Shaking speed: 5 -220 rpm, in steps of 5
 Stroke: 20 mm
 Runtime: 0 - 99 h 59 min / continuous
 Electrical supply: 100 - 240 V, 50 / 60 Hz
 Enclosure protection: IP 21
 Power consumption: max. 4 W
 Heat emission: approx. 0.5 - 2 W
 Ambient temperature: 4°C to 45°C
 Relative humidity: ~ 95 % (non-condensing)
 Dimensions (w x d x h): 250 x 415 x 145 mm
 Weight: 8.1 kg

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Loading Capacity (Examples)

Erlenmeyer flasks 100 ml 9 pieces
 Erlenmeyer flasks 500 ml 2 pieces
 Further loading possibilities see page 47/48

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Universal tray KM Mini
No. 0052 119

Coated tray with drillings for individual loading with spring clamps or test tube racks.
 Shaking platform 200 x 295 mm.

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Universal tray KM
No. 0052 030

Coated tray with drillings for individual loading with spring clamps or test tube racks.
 Shaking platform 356 x 300 mm.

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Description

Well suited for operation in refrigerators or incubators

Basic equipment

Basic device incl. basic platform and rubber mat, mains / charger unit

Technical Data

Motion: orbital
 Shaking platform: 200 x 295 mm
 Max. load: 2 kg
 Shaking speed: 5 - 220 rpm, in steps of 5
 Stroke: 20 mm
 Charging time (battery): 10 hours
 Runtime of battery: min. 40 hours
 Runtime with mains cable: 0 - 99 h 59 min / continuous
 Electrical supply battery charger: 100 - 240 V, 50 / 60 Hz
 Internal electrical supply: 1 battery 12 V / 7Ah
 Enclosure protection: IP 21
 Power consumption: max. 4 W
 Heat emission: approx. 0.5 - 2 W
 Ambient temperature: 4°C to 45°C
 Relative humidity: ~ 95 % (non-condensing)
 Dimensions (w x d x h): 250 x 415 x 145 mm
 Weight: 10.5 kg

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Loading Capacity (Examples)

Erlenmeyer flasks 10 ml 20 pieces
 Erlenmeyer flasks 50 ml 12 pieces
 Further loading possibilities see page 47/48

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Description

Space saving tube roller with gentle tilting and rotating motion

Basic equipment

Basic device with 6 removable plastic rolls and spill tray

Technical Data

Motion:	rotating, tilting motion
Shaking platform:	200 x 300 mm
Dimensions of rolls:	300 x Ø 32 mm
Number of rolls:	6 pc
Max. load:	2 kg
Shaking speed:	5 - 80 rpm, in steps of 5
Angle of tilt:	3° fixed
Stroke at end of rolls:	15 mm
Runtime:	0 - 99 h 59 min / continuous
Electrical supply:	100 - 240 V, 50 / 60 Hz
Enclosure protection:	IP 21
Power consumption:	max. 6 W
Heat emission:	approx. 0.5 - 2 W
Ambient temperature:	5°C to 45°C
Relative humidity:	~ 85 % (non-condensing)
Dimensions (w x d x h):	244 x 465 x 115 mm
Weight:	7.3 kg

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Loading Capacity (Examples)

Suitable for all standard tube sizes, flasks and universal vessels up to 80 mm in diameter thanks to removable rolls

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Overhead Rack


No. 0060 003

with 12 drillings à Ø 29,5 mm (other drillings on request)

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For easy cleaning the plastic rolls (no. 0060 002) as well as the spill tray (no. 0016 416) can be removed.

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Description

Special shaker for small samples in microplates

Basic equipment

Basic device incl. basic platform with holders for 4 microplates

Technical Data

Motion: exact orbital motion
 Shaking platform: 200 x 295 mm for 4 microplates
 Max. load: 2 kg
 Shaking speed: 100 - 1400 rpm, steplessly adjustable
 Stroke: 3 mm
 Timer: 0 - 120 min / continuous
 Electrical supply: 230 V or 115 V, 50 / 60 Hz
 please indicate in case of order
 Enclosure protection: IP 21
 Power consumption: max. 23 W
 Heat emission: approx. 6 W
 Ambient temperature: 5°C to 50°C
 Relative humidity: ~ 85 %
 Dimensions (w x d x h): 220 x 410 x 110 mm
 Weight: 8.5 kg

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Loading Capacity (Examples)

4 microplates
 deepwell-plates, or other plates
 in microplate format

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Shaking platform KM No. 0052 072

For shaker TiMix 2 (instead of standard platform). Flat shaking platform with rim (20 mm) and rubber mat. Shaking platform 200 x 295 mm.

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Slow, three-dimensional motion

Basic equipment

Basic device incl. basic tray and rubber mat

Technical Data

Motion: three-dimensional orbital motion
 Shaking platform: 340 x 340 mm
 Max. load: 10 kg
 Shaking speed: 1 - 23 rpm, steplessly adjustable
 Angle of tilt: 0° to 9° adjustable
 Timer: 0 - 120 min / continuous
 Electrical supply: 230 V or 115 V, 50 / 60 Hz
 please indicate in case of order
 Enclosure protection: IP 21
 Power consumption: max. 18 W
 Heat emission: approx. 6 W
 Ambient temperature: 5°C to 50°C
 Relative humidity: ~ 85 %
 Dimensions (w x d x h): 355 x 450 x 205 mm
 Weight: 18.5 kg

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Multi-storey rack system

No. 0051 247

With 2 trays (340 x 340 mm) and rubber mat.
 For shakers WS 10 and TL 10.
 Height: 200 mm

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(340 x 340 mm)

No. 0051 502

Additional tray for multi-storey rack WS / TL incl. rubber mat.

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Slow, linear tilting motion

Basic equipment

Basic device incl. basic tray and rubber mat

Technical Data

Motion: linear tilting motion
 Shaking platform: 340 x 340 mm
 Max. load: 10 kg
 Shaking speed: 1 - 23 rpm, steplessly adjustable
 Angle of tilt: 5° fixed
 Timer: 0 - 120 min / continuous
 Electrical supply: 115 - 230 V, 50 / 60 Hz
 Enclosure protection: IP 21
 Power consumption: max. 8 W
 Heat emission: approx. 3 W
 Ambient temperature: 5°C to 50°C
 Relative humidity: ~ 85 %
 Dimensions (w x d x h): 355 x 450 x 205 mm
 Weight: 18.5 kg

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The „2 in 1“ shaker

Basic equipment

Basic device incl. rack system with 3 clamping strips h and rubber mat

Technical Data

Motion:	can be switched from orbital to reciprocating (horizontal)
Shaking platform:	300 x 300 mm
Max. load:	2 kg
Shaking speed:	30 - 420 rpm, steplessly adjustable
Stroke:	8 mm
Timer:	0 - 120 min / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 23 W
Heat emission:	approx. 4 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	355 x 455 x 195 mm
Weight:	13.5 kg

Loading Capacity (Examples)

Erlenmeyer flasks 250 ml	9 pieces
Separating funnels 250 ml	1 pieces

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Nr. / No. 6173 000

Description

Reliable, robust shaker for daily use

Basic equipment

Basic device without rack system

Technical Data

Motion KS 15 A:	orbital
Motion KS 15 B:	reciprocating (horizontal)
Shaking platform:	400 x 300 mm
Max. load:	15 kg
Shaking speed:	30 - 420 rpm, steplessly adjustable
Stroke:	17 mm
Timer:	0 - 120 min / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 25 W
Heat emission:	approx. 5 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	510 x 490 x 150 mm
Weight:	18 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	12 pieces
Erlenmeyer flasks 500 ml	6 pieces
Erlenmeyer flasks 1000 ml	4 pieces
Separating funnels 250 ml	1 pieces

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Rack system Combifix KS

No. 0052 071

Rack system with 3 clamping strips h and rubber mat.

Shaking platform 400 x 275 mm, height approx. 80 mm.

For fastening upright vessels. With additional strips also suitable for other vessels.

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Universal tray KS

No. 0051 471

Coated tray with drillings for individual fastening of spring clamps or test tube racks.

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Basic equipment

Basic device without rack system.

Programmable control with touch display

(Technical details of control see page 44)

Technical Data

Motion KS 15 A control:	orbital
Motion KS 15 B control:	reciprocating (horizontal)
Shaking platform:	400 x 300 mm
Max. load:	15 kg
Shaking speed:	30 - 420 rpm, in steps of 5
Stroke:	17 mm
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 25 W
Heat emission:	approx. 7 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	510 x 490 x 150 mm
Weight:	19 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	20 pieces
Erlenmeyer flasks 250 ml	12 pieces
Erlenmeyer flasks 500 ml	6 pieces
Erlenmeyer flasks 1000 ml	4 pieces
Separating funnels 250 ml	1 pieces

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Defined stop*

Nr. 0060 019

Stop always in the same position (KS 15 A control; KS 15 B control on request)

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Nr. / No. 6166 000

Description

High speed for optimal mixing of small sample volumes

Basic equipment


Basic device without rack system

Technical Data

Motion:	exact orbital motion
Shaking platform:	400 x 300 mm for 8 microplates
Max. load:	5 kg
Shaking speed:	100 - 1400 rpm, steplessly adjustable
Stroke:	3 mm
Timer:	0 - 120 min / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 30 W
Heat emission:	approx. 9 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	510 x 490 x 150 mm
Weight:	19 kg

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Standard rack system
No. 0052 101

For 8 standard microplates, microwell or deepwell plates, or other plates of the same size.

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Rack system Combifix KS
No. 0052 071

Rack system with 3 clamping strips h and rubber mat. Shaking platform 400 x 275 mm, height approx. 80 mm. For fastening upright vessels. With additional strips also suitable for other vessels.

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Additional tray
No. 0052 102

As a second storey on the standard rack system TiMix 5. Distance between the 2 trays max. 55 mm.

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Universal tray KS
No. 0051 471

Coated tray with drillings for individual fastening of spring clamps or test tube racks. Shaking platform 430 x 300 mm.

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Rack system with high clamping pins
No. 0052 104

For max. 24 standard microplates, placed directly on top of each other.

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
Rack system with metal clamping pins
No. 0052 125

For max. 48 standard microplates, placed directly on top of each other.

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Nr. 0060 018

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Description

Top technology - also for integration in analytical robot systems

Basic equipment

Basic device without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion:	exact orbital motion
Shaking platform:	400 x 300 mm for 8 microplates
Max. load:	5 kg
Shaking speed:	100 - 1400 rpm, in steps of 10
Stroke:	3 mm (Special stroke on request)
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 30 W
Heat emission:	approx. 11 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	510 x 490 x 150 mm
Weight:	20 kg

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Description

Table-top shaker for high loads and continuous operation

Basic equipment

Basic device with shaking plate and rubber mat, without rack system

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**Rack system
Combifix SM A**
No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels.

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2-storey top frame SM
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 323 mm. Distance between the 2 storeys: 282 mm (height sufficient for 1000 ml Erlenmeyer flasks)

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**Rack system
Combifix SM B**
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v. For fastening horizontal vessels between the clamping strips, or as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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Sliding Plates

Sliding plates with Combifix SM (rack systems Combifix SM, but with additional base plate).

Sliding Plate with Combifix SM A
No. 0051 484

with Combifix SM B
No. 0051 485

with Combifix SM C
No. 0051 486

Mounting set
No. 0052 074

for universal trays, or sliding plates on the 2-storey top frame SM.

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**Rack system
Combifix SM C**
No. 0050 156

Rack system with 2 clamping strips h, 2 clamp strips and 1 spring strip. For fastening separating funnels. Standard version of clamp strip: 5 clamps NS 19 (D20) and 2 clamps NS 29 (D32). With additional strips also suitable for other vessels.

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Universal tray SM
No. 0051 472

Coated tray with drillings for individual fastening of spring clamps or test tube racks. Shaking platform 560 x 400 mm.

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Description

Table-top shaker with programmable shaking tasks for the mixing of difficult media

Basic equipment

Basic device with shaking plate and rubber mat, without rack system.

Programmable control with touch display (Technical details of control see page 44)

Technical Data

Motion SM 30 A control:	orbital
Motion SM 30 B control:	reciprocating (horizontal)
Motion SM 30 C control:	can be switched from orbital to reciprocating
Shaking platform:	560 x 400 mm
Max. load:	30 kg
Shaking speed:	15 - 300 rpm, in steps of 5
Stroke SM 30 A control:	26 mm
Stroke SM 30 B control:	30 mm 46 mm (Option: No. 0052 212) 50 mm (Option: No. 0052 216)
Stroke SM 30 C control:	26 mm
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 46 W
Heat emission:	approx. 8 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	33 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	42 pieces
Erlenmeyer flasks 250 / 300 ml	20 pieces
Separating funnels 250 ml	6 pieces
Separating funnels 1000 ml	2 pieces
Test tube racks	4 pieces

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Defined stop*

No. 0060 020

Stop always in the same position (SM 30 A control and SM 30 B control)

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2-storey top frame SM
No. 0052 065

Top frame with 2 coated trays with drillings for spring clamps or test tube racks. Total height: 323 mm. Distance between the 2 storeys: 282 mm. (height sufficient for 1000 ml Erlenmeyer flasks)

Alternatively, the rack systems Combifix SM or the universal tray SM can be mounted on the 2-storey top frame SM.

For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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Description

Table-top shaker with heating plate and programmable shaking tasks for mixing of difficult media

Basic equipment

Basic device with heating plate, without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion:	orbital
Shaking platform:	560 x 400 mm
Max. load:	25 kg
Shaking speed:	15 - 300 rpm, in steps of 5
Stroke:	26 mm
Runtime:	programmable / continuous
Temperature range:	up to 120°C
Control accuracy:	approx. +/- 1°C
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Enclosure protection:	IP 21
Power consumption:	max. 450 W
Heat emission:	400 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	40 kg

Loading Capacity (Examples)

Erlenmeyer flasks 100 ml	42 pieces
Erlenmeyer flasks 250 / 300 ml	20 pieces
Separating funnels 250 ml	6 pieces
Separating funnels 1000 ml	2 pieces

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*Installation only possible in our works

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**Rack system
Combifix SM B**
No. 0050 155

Basic rack with 4 clamping strips h and 4 clamping strips v as added support for high vessels. The maximum distance between the strips h and v is 60 mm. With additional strips also suitable for other vessels.

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**Rack system
Combifix SM C**
No. 0050 156

Rack system with 2 clamping strips h, 2 cramp strips and 1 spring strip. For fastening separating funnels. Standard version of cramp strip: 5 clamps NS 19 (D20) and 2 clamps NS 29 (D32). With additional strips also suitable for other vessels.

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The below specified shaking speeds against load are approximate values. Depending on specific properties of substances to be shaken these values can differ marginally.

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Description

Programmable shaking tasks for the mixing of difficult media

Basic equipment

Basic device with shaking plate and rubber mat, without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion:	reciprocating (horizontal)
Shaking platform:	560 x 400 mm
Max. load:	30 kg
Shaking speed:	15 - 240 rpm, in steps of 5
Stroke:	50 mm for Gushing Test (No. 0052 216)
Runtime:	programmable / continuous
Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 46 W
Heat emission:	approx. 8 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	680 x 610 x 160 mm
Weight:	33 kg

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2 x Spannleiste h Nr. 0050 400
1 x Spannleiste q Nr. 0050 403

1-storey rack system
No. 0050 154

Basic rack with 5 clamping strips h. For fastening different upright vessels. With additional strips also suitable for other vessels. (see picture)

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2-storey top frame
No. 0052 065

Top frame with 2 coated trays with drillings for mounting rack systems Combifix SM or for spring clamps or test tube racks.

Total height: 323 mm;
Distance between the 2 storeys: 282 mm.

Dimensions (w x d x h):
680 x 420 x 480 mm
Weight: 7.7 kg

Alternatively, the universal tray SM can be mounted on the 2-storey top frame SM. For easy handling (loading away from the shaker) it is recommended to use universal trays SM or sliding plates with Combifix SM, especially on the lower level.

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3-storey top frame
No. 0039 206

Top frame with 3 coated trays for mounting the rack systems Combifix SM. Total height: 400 mm; Distance between the storeys: 177 mm.

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Nr. / No. 6161 000

Description

Constant temperatures; also for high vessels
(e.g. 1000 ml Erlenmeyer flasks)
Can be combined with all compact shakers or used as independent unit

Basic equipment

Coated metal housing with 3 windows, the front door opens upwards,
cable feedthrough Ø 12,6 mm at the back

Technical Data

Temperature range:	5°C above ambient to +50°C
Temperature accuracy:	+/- 1% of the set value
Inner dimensions (w x d x h):	490 x 410 x 310 mm (height sufficient for 1000 ml Erlenmeyer flasks)
Temperature sensor:	PT 100
Air circulation:	240 m ³ /h
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Power consumption:	max. 345 W
External Dimensions (w x d x h):	510 x 450 x 455 mm
Weight:	21 kg

Safety cut-off of the shaker if combined with a Bühler Compact Shaker
Standard temperature controller 5°C above ambient to +50°C

Option

Extended temperature range up to +60°C

No. 0052 105

Installation of an additional heater in order to reach temperatures up
to max. +60°C. Installation only possible in our works.

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Shelf TH 15

No. 0052 126

Additional tray for samples which
have to be temperature controlled
but need no shaking.
Dimensions: 378 x 318 mm.

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Nr. / No. 6162 000

Description

Excellent system thanks to reproducible ambient temperatures and useful accessories
Can be combined with all universal shakers or used as independent unit

Basic equipment

Coated metal housing with 3 windows, the front door opens upwards, socket inside, cable feedthrough Ø 12,6 mm at the back

Technical Data

Temperature range:	5°C above ambient to +50°C
Temperature accuracy:	+/- 1% of the set value
Inner dimensions (w x d x h):	660 x 540 x 430 mm (height sufficient for 2000 ml Erlenmeyer flasks)
Lighting:	standard light 15 W / UV light 15 W
Temperature sensor:	PT 100
Air circulation:	240 m ³ /h
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Power consumption:	max. 620 W
External Dimensions (w x d x h):	680 x 610 x 560 mm
Weight:	33 kg

Safety cut-off of UV light and stop of the shaker when the front door is opened if combined with a Bühler Universal Shaker
Standard temperature controller 5°C above ambient to +50°C

Options

Extended temperature range up to +60°C

No. 0052 105

Installation of an additional heater in order to reach temperatures up to max. +60°C. Installation only possible in our works.

Cooling device TH 30

No. 0051 581

Installation of a cooling coil and connectors for an external cooling system (refrigerated circulator or flow-through cooler). Depending on the cooling capacity of the external cooler temperatures of approx. 5°C below room temperature can be reached (300 - 500 W required). Installation only possible in our works.

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Shelf TH 30
No. 0052 123

Additional tray for samples which have to be temperature controlled but need no shaking.
Dimensions: 365 x 558 mm.

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2-storey top frame SM / TH
No. 0052 117

2-storey top frame for all SM shakers. Dimensions adapted to TH 30. For vessels with a max. height of 14 cm (Erlenmeyer flasks 200 / 250 ml).
Dimensions (w x d x h): 600 x 402 x 192 mm.
For fastening the vessels, either universal trays SM or sliding plates with Combifix SM can be used.

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Description

Big heavy-duty shaker for high loads also during continuous operation
(floor installation only)

Basic equipment

Basic device without shaking plate, without rack system.
Programmable control with touch display
(Technical details of control see page 44)

Technical Data

Motion VKS 75 A control:	orbital
Motion VKS 75 B control:	reciprocating (horizontal)
Shaking platform:	760 x 600 mm
Max. load:	75 kg
Shaking speed:	30 - 200 rpm, in steps of 5
Stroke VKS 75 A control:	26 mm
Stroke VKS 75 B control:	50 mm
	36 mm (Option: No. 0052 213)
	60 mm (Option: No. 0052 214)
	80 mm (Option: No. 0052 215)
	programmable / continuous
Runtime:	
Electrical supply:	100 - 240 V, 50 - 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 70 W
Heat emission:	approx. 15 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	1050 x 835 x 250 mm
Weight:	110 kg

Loading Capacity (Examples)

Erlenmeyer flasks 500 ml	32 pieces
Erlenmeyer flasks 2000 ml	12 pieces
Separating funnels 1000 ml	6 pieces

Further loading possibilities see page 47/48

Options

USB 2.0 interface

No. 0052 931

for the communication between PC and shaker

Control interface

No. 0052 932

for the integration of the shaker in robotic systems (SPS)

Safety Package*

No. 0060 025

Monitoring of motor and shaker speed

RS 232 interface

No. 0052 933

for the communication between robotic system and shaker

Defined stop*

Nr. 0060 021

Stop always in the same position (VKS 75 A control; VKS 75 B control on request)

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**Shaking plate VKS
(1-storey)**

No. 0052 070

Base plate 760 x 600 mm incl. PE plate. Required for mounting rack systems Combifix VKS or universal tray VKS.

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**Universal tray
VKS (1-storey)**

No. 0051 474

Coated tray with drillings for spring clamps or test tube racks. Shaking platform 760 x 600 mm. To be mounted on the shaking plate.

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**Rack system
Combifix VKS A**

No. 0051 487

Basic rack system with 9 clamping strips h. For secure fastening of different upright vessels. With additional strips also suitable for other vessels. To be mounted on the shaking plate.

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**Rack system
Combifix VKS B**

No. 0051 488

Basic rack system with 4 clamping strips h and 4 clamping strips v. For fastening horizontal vessels between the clamping strips, or as added support for high vessels. The maximum distance between the clamping strips h and v is 60 mm. To be mounted on the shaking plate.

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**Rack system
Combifix VKS C**

No. 0051 489

Basic rack system with 2 clamping strips h, 2 clamp strips and 1 spring strip. For fastening separating funnels. Standard version of clamp strip: 7 clamps NS 29 (D32). With additional strips also suitable for other vessels. To be mounted on the shaking plate.

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Total height: 510 mm
Distance
between levels: 387 mm
Weight: 36 kg

Height sufficient for
5000 ml Erlenmeyer flasks

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Total height: 910 mm
Distance
between levels: 340 mm
Weight: 43 kg

Height sufficient for
3000 ml Erlenmeyer flasks

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Platforms for 2- and 3-storey rack systems VKS

Sliding Plate

with Combifix VKS A

No. 0051 493

Base plate 760 x 600 mm with basic
rack and 9 clamping strips h.

with Combifix VKS B

No. 0051 494

Base plate 760 x 600 mm with basic
rack, 4 clamping strips h and 4 clam-
ping strips v.

with Combifix VKS C

No. 0051 495

Base plate 760 x 600 mm with basic
rack and 2 clamping strips h,
2 cramp strips and 1 spring strip.

Universal tray VKS

No. 0051 496

Coated tray with drillings for spring
clamps or test tube racks.
Shaking platform 760 x 600 mm.

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Alternative platforms for 2- and 3-storey rack systems VKS

Sliding Plate

with Combifix SM A

No. 0051 484

Base plate 560 x 400 mm with basic
rack and 5 clamping strips h.

with Combifix SM B

No. 0051 485

Base plate 560 x 400 mm with basic
rack, 4 clamping strips h and 4 clam-
ping strips v.

with Combifix SM C

No. 0051 486

Base plate 560 x 400 mm with basic
rack and 2 clamping strips h,
2 cramp strips and 1 spring strip.

Universal tray SM

No. 0051 472

Coated tray with drillings for spring
clamps or test tube racks.
Shaking platform 560 x 400 mm.

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Basic equipment

Basic device without shaking plate, without rack system.
 Programmable control with touch display
 (Technical details of control see page 44)

Technical Data

Motion:	reciprocating (horizontal)
Shaking platform:	760 x 600 mm
Max. load:	75 kg
Shaking speed:	up to 80 rpm, in steps of 5
Stroke optional:	120 - 150 mm in 10 mm steps please indicate in case of order
Runtime:	programmable / continuous
Electrical supply:	100 - 240 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 90 W
Heat emission:	approx. 15 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	1151 x 835 x 250 mm
Weight:	160 kg

Loading Capacity

Horizontal or vertical barrels up to 80 l capacity

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**Shaking plate VKS
(1-storey)
No. 0052 070**

Base plate 760 x 600 mm incl. PE plate. Required for mounting the rack system no. 0039 001.

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Außenmaße (B x T x H):
780 x 565 x 450 mm
Gewicht: 12 kg

**Special Rack System for barrels
No. 0039 001**

Basic rack with 3 fixing elements. For fastening horizontal or vertical barrels up to 80 l capacity.

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Außenmaße (B x T x H):
800 x 600 x 200 mm
Gewicht: 15 kg

**Special Rack System
for horizontal barrels
No. 0039 002**

Rack system with base plate. For fastening horizontal barrels up to 50 l capacity

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
Technical Data

Motion:	reciprocating (horizontal)
Shaking platform:	760 x 600 mm
Max. load:	75 kg
Shaking speed:	up to 80 rpm, in steps of 5
Stroke:	130 mm
Runtime:	programmable / continuous
Electrical supply:	100 - 240 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Power consumption:	max. 90 W
Heat emission:	approx. 15 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions shaker (w x d x h):	1050 x 835 x 250 mm
Weight:	110 kg
Total weight (Shaker + rack system):	135 kg

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750 x 650 x 370 mm

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Gewicht: 25 kg

**Rack System for
champagne bottles**
No.0052 287

Rack system for fastening up to 36 horizontal champagne bottles.

To show the fixation of the champagne bottles on the rack system, side walls and lid are omitted on the illustration. The rack system is completely closed.

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- Intuitive operation via capacitive touch display for easy operation also with gloves
- Non-reflecting LCD colour display 110 x 67 mm with resolution 800 x 480 dot
- Bright LED background lighting (400cd/m²)
- Display of date and real time
- Sensor for the detection of an incubator hood on the shaker
- Manual operation or automatic operation
- Pause function for the shaking tasks
- Logfile function for the recording of interruptions of the shaking process
- Possibility to change the rotational direction for optimal mixing (orbital motion)
- Automatic operation:
 - 10 programs with 10 segments each, freely definable
 - Optical and acoustic signal at the end of a program
 - Possibility to program a defined time for starting the shaking action
 - Programmable restart after power failure
 - Alpha-numeric program names

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Options

USB 2.0 interface (socket type B) for the communication between PC and shaker No. 0052 931
 Plug-and-play connection with standard terminal programs. Complete remote control of the shaker.

Interface RS 232 for the communication between robotic system and shaker No. 0052 933
 Serial interface for industrial and robotic systems. Shaker completely remote controlled.

Control interface for the integration of the shaker in robotic systems (PLC) No. 0052 932
 - Relais output: shaker is running / stops
 - Relais output: signal alarm / trouble-free operation
 - Relais output: ready for operation
 - Optocoupler input for external start/stop
 - Analog input for setpoint of the speed
 - Analog output for actual speed value

Defined Stop*
 Predefined stop always in the same position. Especially suited for the use of the device in robotic systems.
 (TiMix 5 control, KS 15 A control, SM 30 A control, SM 30 B control, VKS 75 A control; KS 15 B control and VKS 75 B control on request)


Safety Package*
 The speed of both motor and shaker is continuously monitored. In case of deviations, e.g. due to a defective drive belt, a warning message is displayed and an acoustic signal is audible.

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Clamping strip h

Holder with sponge rubber cover for fastening upright vessels.

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Clamping strip v

Holder with sponge rubber cover in combination with clamping strip h.

For fastening horizontal or high vessels.

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Cramp strip

For fastening horizontal separating funnels.

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Spring strip

For securing the stoppers of separating funnels or other horizontal vessels.

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For universal trays. The sizes are related to Erlenmeyer flasks, but are also suitable for round-bottom flasks, beakers, etc.

10 ml	25 ml	50 ml	100 ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml
0009 642	0009 643	0009 644	0009 645	0009 646	0009 647	0009 648	0009 649	0009 653	0009 652

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To be fastened either with clamping strips on rack systems or with a hinged foot on universal trays. Dimensions 265 x 100 mm each.

With the hinged foot it is possible to vary the angle of inclination of the test tubes.

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Abmessung 289 x 100 mm.

Required for universal trays. Size fitting for all test tube racks. Dimensions 289 x 100 mm.

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¹⁾ with 5 clamping strips h (0052 115) additional ²⁾ with 4 clamping strips h (0052 115) additional ³⁾ with 3 clamping strips h (0052 115) additional
⁴⁾ with 2 clamping strips h (0052 115) additional ⁵⁾ with 2 clamping strips h (0050 118) additional ⁶⁾ with 1 clamping strip h (0050 118) additional
⁷⁾ with 3 clamping strips h (0050 118) additional ⁸⁾ with 3 clamping strips h (0050 400) additional ⁹⁾ with 1 clamping strip h (0050 400) additional
¹⁰⁾ with 5 clamping strips h (0050 387) additional ¹¹⁾ with 1 clamping strip h (0050 387) additional

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¹²⁾ with 1 clamp strip (0050 206) and 1 spring strip (0050 207) additional ¹³⁾ Number of pieces depends on shape and dimensions of the separating funnels
¹⁴⁾ with 2 additional clamps NS29 (0001 138) for replacement

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Nr. / No. 6060 000

Description

Compact laboratory heater with fast heating-up rate

Basic equipment

Basic device with support cross

Technical Data

Max. temperature:	700°C
Heating power:	250 W
Heating area:	60 x 60 mm
Max. ambient temperature:	30°C
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Dimensions (w x d x h):	100 x 105 x 105 mm
Weight:	0.7 kg

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Power control LR 1 No. 6070 000

Stepless control of the current
consumption between 0 - 100 %

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Support rod for stand connection

No. 0012 164

Stainless steel rod for integration
of the heater into a support stand,
length approx. 13 cm

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Description

Laboratory heater with power control for larger volumes

Basic equipment

Basic device incl. integrated power control, support cross

Technical Data

Max. temperature:	900°C
Heating power:	800 W
Heating area:	100 x 100 mm
Power control:	stepless control between 0 - 100 %
Max. ambient temperature:	30°C
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Dimensions (w x d x h):	150 x 150 x 170 mm
Weight:	2.5 kg

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Support rod for stand connection

No. 0012 164

Stainless steel rod for integration of the heater into a support stand, length approx. 13 cm

IRB 6 mit Stativgestell, PVC-Schiene und Halteklammern
IRB 6 with support base, horizontal PVC rod and clamps for PVC rod

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Nr. / No. 9066 000

Description

Modular laboratory heater with 6 x IRB 1 in a housing with 6 power control units
Each heater can be individually used and controlled


Basic equipment

6 x IRB 1 with support crosses, 1 x power control unit LR 6

Technical Data

Max. temperature:	each 700°C
Heating power:	6 x 250 W
Heating area:	6 x 60 x 60 mm
Power control:	6 stepless controls between 0 - 100 %
Max. ambient temperature:	30°C
Electrical supply:	230 V, 50 / 60 Hz 115 V version on request
Heating power:	max. 1200 W
Dimensions (w x d x h):	697 x 305 x 122 mm
Weight:	9 kg

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Support base No. 0052 060

Base with 6 support rod (stainless steel), length of the rods approx. 75 cm (without PVC rod and clamps for vessels)

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Horizontal PVC rod for support base No. 0052 091

For fastening clamps for glass vessels

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Examples of customized solutions

With our comprehensive laboratory device program we can offer coherent answers for most demands which arise from practical laboratory work. We believe that our job is even more. A special challenge are the very particular areas of use for which the standard is not enough. Here we design and produce customized solutions. Customer-oriented, quality conscious and tailor made. Some examples of our customized products:

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