IKA

designed for scientists

Dispergierwerkzeuge / Dispersion tools / Outils dispersants

S 25 EC - C - 18 G

S 25 EC - C - 25 G

S 25 EC - C - 18 G - ST

S 25 EC - C - 25 G - ST

S 25 EC - C - 25 F


S 25 EC - T - C - 18 G

S 25 EC - T - C - 25 G


S 25 EC - T - C - 18 G - ST

S 25 EC - T - C - 25 G - ST


S 25 EC - T - C - 25 F


S 25 EC - C ...


S 25 EC - T - C ...


18 G


18 G - ST

25 G

25 G - ST

25 F

Technische Information

Ursprungssprache

DE

4

Technical information

EN


8

Informations Technique

FR

12

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Content

	Page
Explication of warning symbols	8
Safety instructions	8
Product information	9
Application instructions	9
Unpacking	9
Spare parts list	9
Maintenance and cleaning	10
Warranty	10
Technical Data	11

Explication of warning symbols

 DANGER	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 WARNING	Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
 CAUTION	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 NOTICE	Indicates practices which, if not avoided, can result in equipment damage.
 CAUTION	Indicates hazards arising from a hot surface.

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

 **WARNING** Don't touch rotating parts during operation.

 **CAUTION** When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot. Cool the dispersion tool before remove it from the drive.

 **CAUTION** Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There is a risk of splashing and evaporation of liquids; body parts, hair, clothing and jewelry getting caught.

 **NOTICE** Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed.

- The work has to be interrupted immediately if you notice unusual noise and/or increasing emission of liquid from the upper, sidewise rinsing drilling.
- Observe the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Do not use any damaged dispersion tools.
- The dispersion tool is not suitable for continuous operation. Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute. Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispersing causes the medium to heat up.
- Note the minimum/maximum immersion depth in the "Technical data".

- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator (integrated in the shaft tube). Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Product information

Nomenclature:

S 25 EC: (EC = easy clean) suitable exclusively for the **IKA** dispersing drives **T 25 easy clean digital/T 25 easy clean control** and **T 25 mini control**.

T: with temperature measurement

C: ceramic sliding bearing

18/25: Stator diameter

G: Coarse (achievable degree of fineness)

F: Fine (achievable degree of fineness)

ST: saw teeth (on the stator)

- The temperature measurement with the dispersing tool **S 25 EC - T ...** works exclusively with the **IKA** dispersing drive **T 25 easy clean control / T 25 mini control**.
- The dispersing tool **S 25 EC - T...** has an operating hours counter. This works exclusively with the **IKA** dispersing drive **T 25 easy clean control / T 25 mini control**. This drive can detect the operating hours of the dispersing tool and indicates when the next maintenance is required for the dispersing tool.

Note: For more information on temperature measurement and maintenance intervals, please refer to the operating instructions of the **IKA** dispersing drive **T 25 easy clean control / T 25 mini control**.

- All materials are FDA (Food and Drug Administration) conform.

Application instructions

Application: Wet crushing, dispersions, emulsifying (Batch-operation).

Unpacking

Unpacking:

- Please unpack the device carefully.
- Any damage should be notified immediately to the shipping agent (post office, railway network or logistics company).

Delivery scope:

- Dispersion tool according to you order
- Shaft key (see Fig. 1)
- Flat key (see Fig. 1)
- Flat key **25 F** (only with **S 25 EC-(T)-C-25F**, see Fig. 1)
- Mesh bag (for collecting and packaging small parts of the dispersion tool for cleaning, e.g. washing the rotator, clamp in dishwasher.)
- Hook (for hook and hang the service card, see Fig. 2)
- Technical information
- Short instruction.

Spare parts list

Item	Designation
1	Shaft tube
2	Shaft
3	Clamp
5	Ceramic sleeve
7	Rotor
8	Stator
25	Protective cap

Spare parts diagram see front page.

For orders of spare parts, please specify the type of dispersion tool and the designation of the spare part.

Maintenance and cleaning

Cleaning:

For cleaning, the **TS 25-3** cleaning tool (see Fig. 3) can be used. For this purpose, the appropriate cleaning agent (e.g. solvent) for pre-purification and in another, e.g. Surfactants and into the third demineralized water. Then, operate the tool at low speed in the suitable cleaning agents. The remnants of the substance are separated by the flow during operation. Since the tool can subsequently be removed from the drive and brought into the cleaning position (cleaning status, see Fig. 4). Now, the shaft tube of the dispersion tool can also be completely rinsed inside.

If this is not enough, the tool can also be completely disassembled and all parts can then be cleaned separately.

The dispersing tool can be disassembled as described below (item see front side):

- Remove the clamp (3) from the shaft tube (for this purpose, you can insert the small tongue of the flat key into the gap between the two half-shells and then twist the flat key.
- The rotor (7) is unscrewed from the shaft (2) by means of the flat key and shaft key included in the tool set. Use the flat key to hold the rotor and unscrew it by turning the shaft (2) (counterclockwise) with the shaft key.
- Unscrew the stator (8) from the shaft tube (1) (left-hand thread). Insert shaft key into the front-end slot of the shaft tube. By using the flat key the stator can now be screwed off in a clockwise direction.
- Pull the shaft (2) upwards out of the shaft tube (1).

Note: pay attention to the ceramic sleeve (5)! It can fall down when disassembling the rotor (7) and the shaft (2).

The assembly of the dispersion tool in reverse sequence. Pay attention to the ceramic sleeve/bushing (see Fig. 5)!

NOTICE The rotor must not be tightened too much, otherwise the thread will be destroyed.

For cleaning, the **HS 25 EC-DW** cleaning tool (see Fig. 6) can also be used.

Spare parts order:

Spare parts diagram see front page.

For orders of spare parts please specify the serial number on the shaft tube, the type of dispersion tool and the designation of the spare part.

Repair:

Please send the dispersion tool for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the "Decontamination Certificate" form **IKA** or use the download printout of it from **IKA** website: www.ika.com.

If you require servicing, return the instrument in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Technical Data

		S 25 EC - C - 18 G	S 25 EC - C - 18 G - ST	S 25 EC - C - 25 G	S 25 EC - C - 25 G - ST	S 25 EC - C - 25 F
Working range	ml	10 ... 1500		50 ... 2000		50 ... 2000
Stator/Rotor (Ø)	mm	18/12.7	18/13.4	25/17	25/20	25/18
Gap between rotor and stator	mm	0.3	0.25	0.5	0.5	0.5
Max. allowable speed	rpm	25000				
Max. circumferential speed	m/s	16.6	17.5	22.2	26.2	23.6
Min./max. immersion depth	mm	35/150				
Material in contact with medium		stainless steel AISI 316L, ceramic SiC				
pH range		2 ... 13				
Suitable for solvents		yes				
Suitable for abrasive substances		yes				
Max. working temperature	°C	180				
Sterilization methods ①		all methods				
Ultimate fineness, suspensions	µm	10 ... 50		15 ... 50		5 ... 25
Ultimate fineness, emulsions	µm	1 ... 10		1 ... 10		1 ... 5

		S 25 EC - T - C - 18 G	S 25 EC - T - C - 18 G - ST	S 25 EC - T - C - 25 G	S 25 EC - T - C - 25 G - ST	S 25 EC - T - C - 25 F
Working range	ml	10 ... 1500		50 ... 2000		50 ... 2000
Stator/Rotor (Ø)	mm	18/12.7	18/13.4	25/17	25/20	25/18
Gap between rotor and stator	mm	0.3	0.25	0.5	0.5	0.5
Max. allowable speed	rpm	25000				
Max. circumferential speed	m/s	16.6	17.5	22.2	26.2	23.6
Min./max. immersion depth	mm	45/150				
Material in contact with medium		stainless steel AISI 316L, ceramic SiC				
pH range		2 ... 13				
Suitable for solvents		yes				
Suitable for abrasive substances		yes				
Max. working temperature	°C	120				
Temperature measurement		yes				
Temperature measurement accuracy	K	± 0.5				
Sterilization methods ①		all methods except hot air				
Ultimate fineness, suspensions	µm	10 ... 50		15 ... 50		5 ... 25
Ultimate fineness, emulsions	µm	1 ... 10		1 ... 10		1 ... 5

① Sterilization methods:

Allowed procedures	Sterilization techniques
Damp heat	Autoclave to 121°C at 2 bar positive pressure.
Chemical procedures	By germ - killing dissolvent alcohol, phenol, formalin..., disinfectant remainders are to be removed with germ - free water.
Hot air	Germ - killing by hot-air with 160 to 190 °C (approx. 30 min.).

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G: grossier (degré de finesse atteignable)

F: fin (achievable degree of fineness)

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IKA India Private Limited
IKA England LTD.
IKA Works (Thailand) Co. Ltd.
IKA Works Kenya Ltd.