➤ **USAGE, CARE
AND CLEANING**

GRINDING PARTS AND ACCESSORIES!

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CONTENTS

Thank you

for choosing grinding parts or accessories from FRITSCH. As an internationally leading manufacturer of application-oriented laboratory instruments for sample preparation and particle size analysis, we offer you outstanding quality.

Our packaging is designed for world-wide transport; nevertheless, transport damage may occur in individual cases. Please check the product upon delivery or no later than during unpacking.

Grinding parts are wearing parts and are therefore not subject to the normal guarantee. We hope you will understand that we cannot accept complaints regarding grinding parts that have already been used.

Grinding Parts and Accessories

- Ceramic and Agate 3
- Hardmetal Tungsten Carbide..... 4
- Coated 4
- Steel and Tool Steel 5
- Stainless Steel 5
- Tempered Steel 6

Cleaning of Grinding Parts and Accessories

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Grinding Parts and Accessories

CERAMIC AND AGATE

Ceramic and Agate Please handle grinding parts and accessories of agate or ceramic with particular care. In the case of grinding bowls, it is critical that you comply with our recommended grinding set filling, maximum speeds and minimum and maximum grinding material. A sufficient quantity of grinding material is essential to avoid excess wear.

Agate is a natural product and may exhibit varying material qualities and properties depending on the country of origin and storage location. We use exclusively high-quality agate from Brazil. Various colour nuances and inclusions are part of the natural formation process and do not represent quality flaws.

Never cool agate or ceramic grinding parts with liquid nitrogen or dry ice. Avoid rapid temperature changes in heating and cooling. Never heat agate above 100°C. The resulting tension cracks could destroy the grinding parts.

Please ensure that seal surfaces and support surfaces are clean. When using agate and ceramic grinding bowls in the mill, make sure that the grinding bowl components are set 100% level on seal surfaces and in the grinding bowl mount. Incorrect handling can result in tension cracks that will destroy the grinding parts.

Grinding Parts and Accessories

- **HARDMETAL TUNGSTEN CARBIDE**
- **COATED**

Hardmetal Tungsten Carbide Grinding parts and accessories of hardmetal tungsten carbide are very hard and brittle. When using tungsten carbide rotors, for instance, metallic contaminants in the grinding material can break the rotor blades. This is not a quality flaw.

Hardmetal tungsten carbide always contains a certain quantity of cobalt. Cobalt can corrode, resulting in a flecked appearance on the hardmetal surface. This is also not a quality flaw.

To protect the grinding parts, they are coated with an universal corrosion protection product prior to shipping. Please remove this before the first use with typical household cleaning agents or alcohol (e.g. isopropyl alcohol).

Coated Some grinding parts and accessories are available with a hard coating (TiN, WC). The hard coating is only about 4 µm thick and wears off under abrasive wear.

Grinding Parts and Accessories

➤ STEEL AND TOOL STEEL

➤ STAINLESS STEEL

Steel and Tool Steel For assembled grinding parts (e.g. rotor body and blades for cutting rotors) and accessories of steel, corrosion may occur under certain conditions (e.g. working in the alkaline or acidic pH range, working with chemical salts of all kinds, contact with other steels or liquids, etc.) despite high material quality. This also applies to high-alloy steels and is not a quality flaw. Please follow our cleaning instructions.

Our grinding parts and accessories of chromium-free steel, manganese steel and tool steel consist of non-corrosion-resistant steel.

Stainless Steel Grinding parts and accessories of stainless steel have a high chromium content. The precious metal combines with oxygen in the air and creates a very thin, self-regenerating passivation layer on the material surface. However, stainless steel is in no way guaranteed to be rust-free, rather it is rust-resistant steel that offers excellent corrosion protection when cared for properly. Therefore, should corrosion occur, this is not a quality flaw.

To avoid contamination, our corrosion-resistant grinding parts are not treated with additional corrosion protection agents.

Grinding Parts and Accessories

TEMPERED STEEL

Tempered Steel Tempered steels generally have a lower chromium content, meaning that no passivation layer forms. They are characterised by greater hardness rather than corrosion resistance.

To protect the grinding parts, they are coated with an universal corrosion protection product prior to shipping. Please remove this before the first use with typical household cleaning agents or alcohol (e.g. isopropyl alcohol).

Cleaning of Grinding Parts and Accessories

STEEL

Steel Grinding parts may be cleaned with typical metal cleaning agents, whereby mechanical brushing is recommended to support the cleaning process. After cleaning, steel grinding parts and accessories should be dried quickly. Take care to ensure that stainless steels are not stored together with impure steels. To clean such materials, please only use tools that are reserved exclusively for working with stainless steel. If the grinding parts are not used for a prolonged period, we recommend applying corrosion protection oil. Please remove any corrosion flecks that might occur using fleece and steel cleaner, household scouring agents or grinding paste.

Grinding

SAFETY INFORMATION

Safety Information When grinding water suspensions in steel or hard-metal bowls with steel and hardmetal balls, high pressures can arise very quickly. Unforeseen reactions may occur in the grinding bowl. Please always adhere to the safety instructions in the instrument operating manuals.

During grinding, the grinding balls run along the inner wall of the grinding bowl and against the lid. Depending on the size of the grinding balls, grooves may form at the preferred positions, which may vary in severity depending on the load. This is not a quality flaw.

We also refer you to the special safety and operating instructions of the individual mill operating manuals. Our products are subject to constant goods receipt inspections and final inspections in accordance with our ISO 9001 certification.

Please also take note of our element analysis on the accessory pages of our website, www.fritsch.de

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Fritsch GmbH