

SC-96A, block

DESCRIPTION

96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8 - 0.2ml strips or 96 tubes of 0.2 ml.

SPECIFICATIONS

96-well unskirted or semi-skirted microplate (0.2 ml)	1
Max. RPM	1400
Max. temperature	100°C



CAT. NUMBER

BS-010120-FK

SC-96A



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)