




Loading Capacity

CertoClav Vacuum Pro 22L


Tray Set

Order number: 8500481



Bottle Plate

Order number: 8500481A



Article	Number	Width [mm]	Height [mm]	Length [mm]	Quantity
Tray Set	8500481	170	50	430	2
		220	50	430	1
Bottle Plate	8500481A	160	2	400	1

Article	Number	Volume [ml]	Height [mm]	Diameter [mm]	Capacity VAC PRO 22L	Capacity VAC PRO 18L	Capacity VAC PRO 12L	Capacity VAC PRO 8L
Duran Bottle	21801175	50	91	46	24	18	16	12
Duran Bottle	218812854	125	124	55	17	14	11	9
Duran Bottle	218813653	250	158	66	10	8	7	0
Duran Bottle	218814452	500	193	78	8	6	0	0



Certoclav Sterilizer GmbH
Georg Grinninger-Str. 37
A-4050 Traun



Tel: +437229689
E-Mail: support@certoclav.com
Website: www.certoclav.com




Bestellungen, Händlerverifizierung
Sonderanfertigungen
Bedienungsanleitungen und Bedienungsvideos
Technischer Support



Loading Capacity CertoClav Vacuum Pro 18L

Tray Set

Order number: 8500480



Bottle Plate

Order number: 8500480A



Article	Number	Width [mm]	Height [mm]	Length [mm]	Quantity
Tray Set	8500480	170	50	330	2
		220	50	330	1
Bottle Plate	8500480A	160	2	300	1

Article	Number	Volume [ml]	Height [mm]	Diameter [mm]	Capacity VAC PRO 22L	Capacity VAC PRO 18L	Capacity VAC PRO 12L	Capacity VAC PRO 8L
Duran Bottle	21801175	50	91	46	24	18	16	12
Duran Bottle	218812854	125	124	55	17	14	11	9
Duran Bottle	218813653	250	158	66	10	8	7	0
Duran Bottle	218814452	500	193	78	8	6	0	0




Certoclav Sterilizer GmbH
Georg Grinninger-Str. 37
A-4050 Traun



Tel: +437229689
E-Mail: support@certoclav.com
Website: www.certoclav.com




Bestellungen, Händlerverifizierung
Sonderanfertigungen
Bedienungsanleitungen und Bedienungsvideos
Technischer Support



Loading Capacity CertoClav Vacuum Pro 12L

Tray Set

Order number: 8500479



Bottle Plate

Order number: 8500479A



Article	Number	Width [mm]	Height [mm]	Length [mm]	Quantity
Tray Set	8500479	140	40	340	2
		180	40	340	1
Bottle Plate	8500479A	125	2	320	1

Article	Number	Volume [ml]	Height [mm]	Diameter [mm]	Capacity VAC PRO 22L	Capacity VAC PRO 18L	Capacity VAC PRO 12L	Capacity VAC PRO 8L
Duran Bottle	21801175	50	91	46	24	18	16	12
Duran Bottle	218812854	125	124	55	17	14	11	9
Duran Bottle	218813653	250	158	66	10	8	7	0
Duran Bottle	218814452	500	193	78	8	6	0	0



Certoclav Sterilizer GmbH
Georg Grinninger-Str. 37
A-4050 Traun



Tel: +437229689
E-Mail: support@certoclav.com
Website: www.certoclav.com



Bestellungen, Händlerverifizierung
Sonderanfertigungen
Bedienungsanleitungen und Bedienungsvideos
Technischer Support




Loading Capacity

CertoClav Vacuum Pro 8L

Tray Set

Order number: 8500478



Bottle Plate

Order number: 8500478A



Article	Number	Width [mm]	Height [mm]	Length [mm]	Quantity
Tray Set	8500478	120	33	300	2
		140	33	300	1
Bottle Plate	8500478A	105	2	280	1

Article	Number	Volume [ml]	Height [mm]	Diameter [mm]	Capacity VAC PRO 22L	Capacity VAC PRO 18L	Capacity VAC PRO 12L	Capacity VAC PRO 8L
Duran Bottle	21801175	50	91	46	24	18	16	12
Duran Bottle	218812854	125	124	55	17	14	11	9
Duran Bottle	218813653	250	158	66	10	8	7	0
Duran Bottle	218814452	500	193	78	8	6	0	0



Certoclav Sterilizer GmbH
Georg Grinninger-Str. 37
A-4050 Traun



Tel: +437229689
E-Mail: support@certoclav.com
Website: www.certoclav.com



Bestellungen, Händlerverifizierung
Sonderanfertigungen
Bedienungsanleitungen und Bedienungsvideos
Technischer Support