

BRANDnew! Transferpette® pro

Available from September 2025



Single-channel micropipette Transferpette® pro, adjustable, DE-M

	Capacity, µl (color-coded)			A* ± %	µl	CV* ≤ %	µl	Subdivision µl	Cat. No.
●	0.1	–	1	2	0.02	1.2	0.012	0.001	706868
●	0.1	–	2.5	1.4	0.035	0.7	0.018	0.002	706869
●	0.5	–	10	1	0.1	0.5	0.05	0.01	706870
●	2	–	20	0.8	0.16	0.4	0.08	0.02	706871
●	2	–	20	0.8	0.16	0.4	0.08	0.02	706872
●	5	–	50	0.8	0.4	0.3	0.15	0.05	706873
●	10	–	100	0.6	0.6	0.2	0.2	0.1	706874
●	20	–	200	0.6	1.2	0.2	0.4	0.2	706878
●	30	–	300	0.6	1.8	0.2	0.6	0.5	706879
●	100	–	1.000	0.6	6	0.2	2	1	706880
●	250	–	2.500	0.6	15	0.2	5	2	706881
●	500	–	5.000	0.6	30	0.2	10	5	706882
●	1.000	–	10.000	0.6	60	0.2	20	10	706884



NEW!

Multi-channel micropipette Transferpette® pro, adjustable, DE-M

	Capacity, µl (color-coded)			A* ± %	µl	CV* ≤ %	µl	Subdivision µl	Cat. No.
●	0.5	–	10	1.6	0.16	1	0.1	0.01	706900
●	5	–	50	0.8	0.4	0.4	0.2	0.05	706906
●	10	–	100	0.8	0.8	0.3	0.3	0.1	706908
●	20	–	200	0.8	1.6	0.3	0.6	0.2	706910
●	30	–	300	0.8	2.4	0.3	0.9	0.5	706912
●	50	–	1.250	1	12.5	0.3	3.75	1	706916



Multi-channel micropipette Transferpette® pro -12, adjustable, DE-M

	Capacity, µl (color-coded)			A* ± %	µl	CV* ≤ %	µl	Subdivision µl	Cat. No.
●	0.5	–	10	1.6	0.16	1	0.1	0.01	706920
●	5	–	50	0.8	0.4	0.4	0.2	0.05	706926
●	10	–	100	0.8	0.8	0.3	0.3	0.1	706928
●	20	–	200	0.8	1.6	0.3	0.6	0.2	706930
●	30	–	300	0.8	2.4	0.3	0.9	0.5	706932
●	50	–	1.250	1	12.5	0.3	3.75	1	706936



* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth operation. A = Accuracy, CV = Coefficient of variation

BRAND. For lab. For life.®