



NEW!



Transferpette® *pro*

Your tasks. Your pipette.
Precise, ergonomic, and flexible.

BRAND. For lab. For life.®

Transferpette® *pro* micropipette

The Transferpette® *pro* micropipette allows you to work precisely, ergonomically, and flexibly – just as your daily lab routine demands!

Its slim design and adjustable finger rest ensure that the pipette fits perfectly in your hand, making it ergonomic and comfortable to use, even during long pipetting sessions.

To keep you flexible at all times, the Transferpette® *pro* features two adjustment techniques: Easy Calibration for permanent calibration without tools and User Adjustment for quick, reversible adaptation to changing conditions such as temperature, tip shape, or media properties.

Therefore: Your tasks. Your pipette.

Your benefits:

- + Simple: Adjustable finger rest and slim design for more comfort
- + Easy: User Adjustment for changing conditions
- + Efficient: Reliable, durable and easy to maintain

Your tasks. Your pipette.

Every application has different requirements. With the Transferpette® *pro*, you are fully prepared for this, as it adapts perfectly to your individual requirements with additional innovative features:

Adjustable finger rest



Adjustable to 9 positions for an ergonomic hand position

User Adjustment



Easy adjustment to special conditions such as temperature, tip shape, or media properties

Counter protection



Reliable protection against over-rotation beyond the mechanical limit

Extended volume range



Now also with 300 µl, 2.5 ml, and 1250 µl nominal volume limit



Counter protection

Longevity through additional protection of the mechanics

Adjustable finger rest

Simple adjustment without tools

User Adjustment

Temporary adjustment to changing conditions

Extended volume range

13 single-channel pipettes from 0.1 µl to 10 ml
6 multi-channel pipettes from 0.5 µl to 1250 µl



All product information can be found at shop.brand.de

Transferpette® *pro* and Transferpette® *S* in comparison

	Transferpette® <i>pro</i>	Transferpette® <i>S</i>
Volume range single-channel (adjustable)	12 sizes: from 0.1 µl to 10 ml	10 sizes: from 0.1 µl to 10 ml
Volume range multi-channel	6 sizes: from 0.5 µl to 1250 µl	5 sizes: from 0.5 µl to 300 µl
Size of handle	Slim	Standard
Adjustable finger rest	✓	-
One-handed operation	✓	✓
Universal tip cone	✓	✓
Volume adjustment protection	✓	✓
Counter protection	✓	-
Easy Calibration Technology	✓	✓
User Adjustment	✓	-
Completely autoclavable (121 °C)	✓	✓
UV-resistant	✓	✓
Warranty for single-channel pipettes	3 years	3 years

