



designed for scientists



C 248 Oxygen filling station

/// Data Sheet

The C 248 oxygen filling station with manometer is designed to fill the C 5010 decomposition vessel with oxygen and is generally included in all C 200 calorimeter packages.

The scope of delivery also includes a flexible two meter pressure hose in order to improve the positioning of the station. In case the C 200 calorimeter or the C 248 oxygen filling station need to be operated with a C 5012 decomposition vessel, the C 200.2 conversion kit is required. This accessory is sold separately.



designed for scientists

Technical Data

Weight [kg]	2.3946
-------------	--------